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ULTIMATE  
BODY  
COMPOSITION  
ANALYZER

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
## InBody 270 Body Composition Analysis Manual

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Remove the socks and step on the footplate

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
Wait for the weight measurement finish then  
enter your phone number (professional mode)  
or age and height (self mode)

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



The results will be shown on screen while the test finished and wait for print out

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個人資料及測試日期		ID	Height	Age	Gender	Test Date & Time
		Jane Doe	156.9cm	51	Female	2018.04.03. 09 : 46
身體成份分析						
Body Composition Analysis						
1	Total amount of water in body	Total Body Water	(L)	26.5 ( 26.4 ~ 32.2 )		
2	For building muscles	Protein	(kg)	7.2 ( 7.0 ~ 8.6 )		
3	For strengthening bones	Minerals	(kg)	2.64 ( 2.44 ~ 2.98 )		
4	For storing excess energy	Body Fat Mass	(kg)	22.8 ( 10.3 ~ 16.5 )		
5	Sum of the above	Weight	(kg)	59.1 ( 43.9 ~ 59.5 )		

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# Before using InBody

Few groups of people are not allowed to use InBody :


- Pregnant women
- People with implanted medical electronic devices, such as pacemakers
- People who have joint replacement surgery
- When a child or people with disabilities is using InBody, must be under supervision or assistance.

Before step on InBody, please make sure to remove all heavy items and metal objects on your body, then you can start the measuring.

\*Assessment is for reference only not for any diagnostic purpose


# InBody 270 Body Composition Analysis Manual

## 1. Please remove the socks and step on the footplate


## 2. The screen will show the weight

Wait for the weight measurement finish then enter your phone number (professional mode) or age and height (self mode)

(Professional mode)

(Self mode)


## 3. Grip and hold the hand electrodes

Maintain proper posture until the test finish

# InBody 270 Body Composition Analysis Manual

## 4. Wait for the measurement finish


(Professional mode)


(Self mode)


## 5. The results will be shown on screen while the test finished and wait for print out

(Professional mode)

(Self mode)

# InBody Report Interpretation

Personal information  
& Test Date

Body Composition  
Analysis

ID	Height	Age	Gender	Test Date & Time
Jane Doe	156.9cm	51	Female	2018.04.03. 09 : 46

Body Composition Analysis			
1	Total amount of water in body	Total Body Water	(L) 26.5 ( 26.4 ~ 32.2 )
2	For building muscles	Protein	(kg) 7.2 ( 7.0 ~ 8.6 )
3	For strengthening bones	Minerals	(kg) 2.64 ( 2.44 ~ 2.98 )
4	For storing excess energy	Body Fat Mass	(kg) 22.8 ( 10.3 ~ 16.5 )
5	Sum of the above	Weight	(kg) 59.1 ( 43.9 ~ 59.5 )

## 1 Total Body Water

- Body water content

## 2 Protein

- One of the component build up your muscle

## 3 Minerals

- Which can maintain body health & metabolism

## 4 Body Fat Mass

- Body fat weight

## 5 Weight

- Sum of all values

Values in parentheses are ideal ranges based on your measurements

# InBody Report Interpretation

## Muscle-Fat Analysis

### 1. Weight

### 2. Skeletal Muscle Mass

Muscle responsible for movement and posture ; most easily affected through exercise.

		Under		Normal		Over						
	(kg)	55	70	85	100	115	130	145	160	175	190	205
Weight	(kg)	66.4										
SMM <small>Skeletal Muscle Mass</small>	(kg)	26.7										
Body Fat Mass	(kg)	18.1										

### 3. Body Fat Mass

The value will divide into three range

Under

Normal

Over

		Under		Normal		Over						
	(kg)	55	70	85	100	115	130	145	160	175	190	205
Weight	(kg)	66.4										
SMM <small>Skeletal Muscle Mass</small>	(kg)	26.7										
Body Fat Mass	(kg)	18.1										

- >> InBody compares your results to other people who have the same [height](#) and [gender](#)
- >> When the value reaches 100%, it reflects your weight is equivalent to a healthy average
- >> If the value reaches 130%, it reflects your weight is **above** the average 30%
- >> In contrast, if the value reaches 70%, it reflects your weight is **below** the average 30%



		Under		Normal		Over						
	(kg)	55	70	85	100	115	130	145	160	175	190	205
Weight	(kg)	66.4										
SMM <small>Skeletal Muscle Mass</small>	(kg)	26.7										
Body Fat Mass	(kg)	18.1										

# InBody Report Interpretation

## Obesity Analysis

### 1. BMI

$$\frac{\text{Weight (kg)}}{\text{Height}^2 (\text{m}^2)}$$


- Underweight: <18.5
- Normal: 18.5-22.5
- Risk to overweight: 22.6-24.9
- Overweight: 25-29.9
- Obese: ≥30

### 2. Percent Body Fat (PBF)

$$\frac{\text{Body fat}}{\text{Mass}} \times 100\%$$

Body weight

Body Composition Analysis			
Total amount of water in body	Total Body Water	(L)	26.5 ( 26.4 ~ 32.2 )
For building muscles	Protein	(kg)	7.2 ( 7.0 ~ 8.6 )
For strengthening bones	Minerals	(kg)	2.64 ( 2.44 ~ 2.98 )
For storing excess energy	Body Fat Mass	(kg)	22.8 ( 10.3 ~ 16.5 )
Sum of the above	Weight	(kg)	59.1 ( 43.9 ~ 59.5 )

In **Male**, the healthy range of PBF is between 10-20%

In **Female**, the healthy range of PBF is between 18-28%

# InBody Report Interpretation


## Segmental Lean Analysis

InBody segmental lean analysis provides 5 body parts results :

- 1
- Right Arm
- 2
- Left Arm
- 3
- Trunk
- 4
- Right Leg
- 5
- Left Leg

Example: Left Arm


1.71 kg lean muscle in left arm  
In Normal Range

Each part lean muscle weight

It determines if you have enough muscle mass to support your weight

It determines if you in the healthy range

# InBody Report Interpretation


## Segmental Lean Analysis

InBody segmental lean analysis provides 5 body parts results :

- 1
- Right Arm
- 2
- Left Arm
- 3
- Trunk
- 4
- Right Leg
- 5
- Left Leg

Example: Right Arm


1.7 kg fat muscle in right arm  
In Over Range

Body fat of each part

Compared with the people who have same height and gender. If you have 192% means 92% over the normal range


It determines if you in the health range

# InBody Report Interpretation

## Body Composition History

Graphically tracks the historical trending of your body compositional changes

including:

- ➔ Weight
- ➔ Skeletal Muscle Mass (SMM)
- ➔ Percent Body Fat (PBF)

# InBody Report Interpretation

## InBody Score

Total score that reflect the evaluation of body composition. 70-90 score means that the health status is between “good” and “very good” state. A muscular person may score over 100

<b>InBody Score</b>
<b>81</b> / 100 Points
<small>* Total score that reflects the evaluation of body composition. A muscular person may score over 100 points.</small>

## Weight Control

Is a recommendation that how to achieve ideal weight after body composition analysis

“+” = increase weight

“-” = decrease weight

Example : need to lose 3.7kg fat and increase 3.5 kg muscle

\*InBody do not recommended decrease muscle weight

<b>Weight Control</b>	
Target Weight	51.7 kg
Weight Control	- 7.4 kg
Fat Control	- 10.9 kg
Muscle Control	+ 3.5 kg

# InBody Report Interpretation


## Obesity Evaluation

**BMI divided into**

Under | Normal | Slightly Over | Over

**PBF divided into**

Normal | Slightly Over | Over

## Waist-Hip Ratio & Visceral Fat Level


### Waist-Hip Ratio = $W \div H$

**Female** ideal waist-hip ratio : 0.75-0.85

**Male** ideal waist-hip ratio : 0.8-0.9

If it exceeds the range, it defined as central obesity

High visceral fat linked with an increased risk of developing cardiovascular disease.

# InBody Report Interpretation

## Research Parameters

### 1. Fat Free Mass

- Body Weight without fat  
Fat Free Mass=  
Body Weight-Body Fat Mass

Research Parameters	
1	Fat Free Mass 36.3 kg
2	Basal Metabolic Rate 1154 kcal (1255~1451)
3	Obesity Degree 114 % ( 90~110 )
4	SMI 5.8 kg/m <sup>2</sup>

### 2. Basal Metabolic Rate

- Minimum amount of energy, or calories, required to sustain life at a resting state within a 24-hours period

### 3. Obesity Degree

- is a percentage of ideal body weight the normal range is 10% above or below (i.e. 90-110%)

### 4. Skeletal Muscle Mass Index, SMI

- Skeletal muscle mass index was calculated by dividing the limb skeletal muscle mass by the square of the height

# InBody Report Interpretation

## Calorie Expenditure of Exercise

ID	Height	Age	Gender	Test Date & Time
Jane Doe	156.9cm	51	Female	2018.04.03. 09 : 46

Body Composition Analysis			
Total amount of water in body	Total Body Water	(L)	26.5 ( 26.4 ~ 32.2 )
For building muscles	Protein	(kg)	7.2 ( 7.0 ~ 8.6 )
For strengthening bones	Minerals	(kg)	2.64 ( 2.44 ~ 2.98 )
For storing excess energy	Body Fat Mass	(kg)	22.8 ( 10.3 ~ 16.5 )
Sum of the above	Weight	(kg)	59.1 ( 43.9 ~ 59.5 )

Calorie Expenditure of Exercise			
Golf	115	Gateball	131
Walking	148	Yoga	148
Badminton	150	Table Tennis	150
Tennis	197	Bicycling	197
Boxing	197	Racketball	197
Hiking, No load	229	Jumping Rope	235
Aerobics	235	Jogging	235
Soccer	235	Swimming	235
Japanese Fencing	250	Racketball	250
Squash	250	Taekwondo	250

↗Based on your current weight  
 ↗Based on 30 minute duration

The calorie expenditure after 30mins exercise,  
base on the participant's weight.

Example : 59.1 kg Female

After jumping rope, aerobics, swimming,  
football or jogging

for 30mins can expend 235 calories

# InBody Q&A

## 1. Why the measurement is not successful?

During the measurement, the hand or foot moves away from the electrode contact point. Also, dry skin will cause measure error. It is recommended to thoroughly wipe the palms and the feet with a towelette before taking the measurement.

## 2. Why does eating affect measurement accuracy?

After a meal, the food and water in the stomach may be considered as fat. The stomach wall interferes with the flow of electricity. Body weight can easily change by 1~2kg due to eating or excreting. Therefore, it is not recommended to test immediately after a meal.

## 3. Why we cannot exercise before the measurement?

Exercise will change the distribution of water in the body. You can feel that the muscles become stiffness, since the water is concentrated in the muscles. The result may not accurately. Therefore, you should take measurements at least 2~3 hours after exercise.

# InBody Q&A

## 4. Why the result is not the same when test again within 5 minutes ?

Participants may not strictly assume the exact same posture for each test. Moreover, during the rest period, the water distribution of the legs is easy to change, and it is not easy to control the difference. Therefore, it is recommended to wait another 1~2 minutes before taking another test. It is not recommended to test immediately after walking or walking stairs.

## 5. How often should I take the InBody Test?

When the weight management program is in progress, it is recommended to take another InBody test after 2 weeks. You should have another InBody test within 4 weeks at the latest.

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