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Before Using PASESA

**PASESA is only suitable for
the people who are over 18 years old**

Few groups of people are not allowed to use PASESA :

- Pregnant women
- People with wound, swelling or bruises on their arm
- People with implanted medical electronic devices, such as pacemakers, electronic life support systems, such as artificial hearts/lungs, portable medical electronic products, such as electrocardiographs are not allowed to use

Before Measurement



First, please check if the devices are properly connected (including the power supply)

If the measurement cannot be started/paused when measurement / the result cannot be displayed after the measurement

1. Please double check if the devices are properly connected (including the power supply)
2. Start measurement again
3. If still can't measure properly, please restart the device and try measurement again
4. If still can't measure properly, please contact our staff

PASESA setting

1. Open the elbow rests and the arm barrel facing you.
2. Put PASESA near the left side of your body.
3. The height difference between tables and chairs is **25-30 cm**. Inappropriate height difference will cause the body oblique.
(The heart should be at the same level as the sign)
4. It is recommended to keep the distance between the machine and the table about 2 cm and the distance between the body and the table about 2 cm.
(Please keep your back straight and have appropriate distance)
5. Inappropriate posture during measurement will cause inaccurate result
(Bend over, sitting cross-legged or low table, etc.)

PASESA Measurement Procedure

1. Rest for a few minutes to keep your mind relaxed. Please take off your heavy clothes and put down your backpack.
2. Sit straight, put your left hand through the arm barrel and palm up. Place your elbow at the groove of the elbow rest.
(The heart should be at the same level as the sign)
3. Please do not talk, laugh, move your body during the measurement.
4. Please keep your body and mind calm during the measurements.
5. The measurement lasts long about 2 minutes. PASESA will have prompt sound when the measurement is finished.
6. Please wait until the cuff is completely released and the staff indicates that the measurement is completed, then withdraw the arm and move the body.
7. Measurement finished.

Feeling unwell during the measurement

- The arm has numbness during the measurement is normal situation.
- If the participant feels discomfort during the measurement, please notify the staff to press the stop button.
- During the measurement, if you feel lack of blood circulation in the fingers. Please raise up your hands above your head and squeeze and stretch your palms about 15 times to relieve this symptoms.


PASESA Report Interpretation

Arterial velocity-pulse index (AVI)
An index of arterial reflected waves. It determines the arteriosclerosis level in the central artery

Central Systolic Blood Pressure (CSBP) refers to the pressure in the aorta near the heart. Lower CSBP helps to prevent cardiovascular disease risk.
Normal Range : 80-130

Central Aortic Pulse Pressure (CAPP) is increasingly considered a better cardiovascular prognostic marker than conventional cuff brachial blood pressure.
Normal Range : 25-49

Arterial pressure-volume index (API)
An index of volume of a conductive blood vessel. It determines the arteriosclerosis level in the peripheral artery

	Systolic Pressure	Diastolic Pressure
Normal	Under 120	Under 80
Elevated	120-139	80-89
High Blood Pressure	Above 140	Above 89

Pulse, a rhythmical throbbing of the arteries as blood is propelled through them.

- ➔ Adult normal pulse range : 60-100/min
- ➔ Systolic blood pressure indicates how much pressure your blood is exerting against your artery walls when the heart beats
- ➔ Diastolic blood pressure indicates how much pressure your blood is exerting against your artery walls while the heart is resting between beats.

Arterial velocity-pulse index (AVI)

An index of arterial reflected waves.

It determines the arteriosclerosis level in the central artery

High risk of central atherosclerosis when AVI value get higher

AVI value easily affected by medicine, Chinese medicine and lifestyle

Age	Lower Limit	Upper Limit	Slightly high ($\leq 20\%$ above the limit)	Moderately high (Between 20%-50%)	Severely high (more than 50%)
20-24	10	16	17	20	24
25-29	11	18	19	22	27
30-34	12	18	19	22	27
35-39	13	20	21	24	30
40-44	14	22	23	27	33
45-49	15	23	25	28	35
50-54	16	25	27	30	38
55-59	17	27	29	33	41
60-64	18	30	32	36	45
65-69	20	32	34	39	48
70-74	22	35	37	42	53
75-79	23	38	40	46	57
80	25	39	41	47	59

Arterial pressure-volume index (API)

An index of volume of a conductive blood vessel.

It determines the arteriosclerosis level in the peripheral artery.

High risk of peripheral atherosclerosis when AVI value get higher

AVI value easily affected by medicine, Chinese medicine and lifestyle

Age	Lower Limit	Upper Limit	Slightly high ($\leq 20\%$ above the limit)	Moderately high (Between 20%-50%)	Severely high (more than 50%)
20-24	15	22	24	27	33
25-29	16	23	25	28	35
30-34	16	24	26	29	36
35-39	16	25	27	30	38
40-44	17	26	28	32	39
45-49	18	27	29	33	38
50-54	18	27	29	33	38
55-59	18	29	31	35	41
60-64	19	30	32	36	45
65-69	20	31	33	38	47
70-74	21	32	34	39	53
75-79	21	33	35	40	59
80	22	33	35	40	59

When AVI and API values are high, the following actions are recommended

Value	Detail explanation	Action
Slightly high ($\leq 20\%$ above the limit)	We recommend that you should have a cardiovascular examination. Also, develop healthy lifestyle, such as regular exercise, low-fat, low-sugar, low-sodium, and high-fiber diet, and control your weight to avoid further aggravation of arteriosclerosis	Having Pasesa review after 3 months
Moderately high (Between 20%-50%)	We recommend that you should have a comprehensive cardiovascular examination, and have treatment according to the doctor's advice	Having Pasesa review and adjust the treatment plan after 3 months
Severely high (more than 50%)	It generally indicates that there have cardiovascular disease. You should have medical treatment as soon as possible	Need to consult doctor as soon as possible

Abnormal AVI & API value explanation

	Normal	Upper Limit	Lower Limit	≤20% above the limit	Above the limit between 20%-50%	Above the limit more than 50%
AVI		High		The value is slightly higher. We recommend that you should develop healthy lifestyle, such as regular exercise, low-fat, low-sugar, low-sodium, and high-fiber diet, and control your weight to avoid further aggravation of arteriosclerosis	We recommend that you should develop a healthy lifestyle. Also, go to the hospital for comprehensive cardiovascular examinations regularly and receive preventive treatment based on the results to reduce the risk of serious cardiovascular diseases	It indicates that the risk of suffering from cardiovascular diseases (such as coronary heart disease, cerebral apoplexy) in the next few years is high, and need to have detailed cardiovascular examinations and appropriate treatment according to the doctor's advice
API		High				

AVI			Low	People who exercise regularly and have a healthy lifestyle. The value may be lower, since the blood vessel elasticity is better than normal people. However, in the perspective of Chinese medicine, it is a symptom that people are a deficiency of "Qi" and blood. They will appear weak and sick, pale and weak pulse, etc. In addition, for some people with low blood pressure, the measured value may be lower than the lower limit.
API			Low	

AVI	Normal			The API value will be affected by the changes in blood pressure or sympathetic nerve activity. Psychological tension will cause high API value It is recommended to take another measurement after rest. In addition, some studies have shown that upper arm resistance training will cause upper arm arterial stiffness. Therefore, people who regularly have arm resistance training may also have higher API value
API		High		

AVI		High		AVI and API reflect the stiffness of central and brachial arteries. The central artery stiffness is more closely related to the risk of cardiovascular disease. If the AVI values higher than the upper limit, it indicates a high risk of cardiovascular disease, and we recommend that you should go to the hospital for cardiovascular examination.
API	Normal			

AVI	Normal			AVI and API reflect the stiffness of central and brachial arteries. The central artery stiffness is more closely related to the risk of cardiovascular disease. If the AVI values higher than the upper limit, it indicates a high risk of cardiovascular disease, and we recommend that you should go to the hospital for cardiovascular examination.
API			Low	

AVI			Low	People who have low blood pressure, the central artery may be softer than normal people. It may cause low AVI value. Also, people who have arrhythmia will cause inaccurate measurement, we recommend that you should take another measurement after rest.
API	Normal			

PASESA Abnormal Result

- Arm circumference over or less than 19 - 32 cm
- Inappropriate posture during measurement
(Bend over, Sitting on a lower chairs, etc.)
- People who are taking medicine or feeling unwell
- Other factors such as measurement time, clothing, mood changes, etc.
- Please do not eat, drink alcohol, smoke, exercise or bath before measurement 30 minutes
- Please empty your bladder before measurement.
- Please take off thicker clothes such as coats and sweaters.
(Please do not measure when you rolling up your sleeves)
- Please rest more than 10 minutes and measure in a quiet and suitable environment

The result cannot appear properly after measurement

01

1. If the customer has arrhythmia, the result will not appear properly.

02

1. First, please doublecheck if the devices are properly connected (including the power supply)
2. Start measurement again
3. If still can't measure properly, please restart the device and try measurement again
4. If still can't measure properly, please contact our staff

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