

Galvanic Spa™ System II and Galvanic Spa System II EX

BRINGS THE SPA TO YOU



Product Overview

Expand your home spa experience. Rejuvenate your complexion, revitalize your scalp, and refine your body with professional spa treatments—no appointment necessary. Now you can bring more of the spa experience home with the patented Galvanic Spa System II.

Small, smart, and powerful, the programmable Galvanic Spa Instrument with patented self-adjusting galvanic currents and interchangeable heads—for the face, scalp, and body—works synergistically with specially formulated products to facilitate the transport of key ingredients for optimal performance. Your radiant complexion; fuller, healthier looking hair; and improved appearance will leave people wondering how you find the time.

Target Audience

Primary:

- Skin-involved, image conscious women ages 25–65 who are interested in spa quality treatments but struggle to justify the time and expense associated with visiting a day spa.

Secondary:

- Men ages 25–65 who want healthy, abundant hair.
- Men ages 25–65 who are image conscious and somewhat skin involved.

Benefits

- Treatments with the Galvanic Spa System Instrument lasting five minutes or more can enhance the delivery of key ingredients for up to 24 hours.
- The Galvanic Spa System features a proprietary, patented instrument with four interchangeable conductors:

Face Conductor—offers all of the benefits of the original fixed face conductor, while adding the flexibility of being removable for other conductor use and other treatment applications.

Scalp Conductor—this comb-like conductor is easy to pull through your hair, maintaining crucial contact between the conductor and your scalp.

Body Conductor—three rounded nodes maximize surface area contact for wider coverage and better massaging action.

Spot Treatment Conductor—this conductor head utilizes a flat, smooth, oval shaped treatment node that centralizes surface

contact area. This node provides better contact with the skin and helps deliver a more concentrated treatment.

Other Features

- **Replaceable batteries**—eliminates cords, adapters, and dependence on electrical outlets.
- **User-selectable polarity**—increases the number of products for use with the Galvanic Spa System.
- **Flexible program timing**—accommodates a wide range of applications: scalp, face, body, and spot treatments.
- **User-friendly handheld contact points**—hand is in a natural, comfortable position.
- **Automatically adjusts to your skin**—delivers the correct amount of current for the most effective treatment.
- **Two-year warranty**—from original date of purchase.

Demonstration

Before performing this magnet demonstration, explain that all galvanic instruments use electrical charges to facilitate the transport of key ingredients for maximum benefit. An important characteristic of electrical charges is utilized by these instruments—like charges repel each other and opposite charges attract each other.

(+) ◀▶ (+) (–) ◀▶ (–)
(+) ▶◀ (–) (–) ▶◀ (+)

If a key ingredient is negatively charged, the negative cycle of the Galvanic Spa Instrument will repel this ingredient—facilitating its transport. In the same way, if the ingredient is positively charged, the instrument will facilitate transport when the cycle is positive. If there are impurities within the skin that need to be pulled out, they will be extracted with an opposite charge.

To demonstrate how Galvanic Spa System currents work synergistically with specially formulated products to facilitate the transport of select key ingredients, use magnets on a magnetic board.

Follow the steps below:

1. Obtain a magnetized board and several magnets (it's best if the magnets are wheel-shaped so they will roll across the board).
2. Draw or paste a picture of a face onto the magnet board.
3. Place one of the wheel-shaped magnets on the edge of the face.

Galvanic Spa™ System II and Galvanic Spa System II EX

4. Connect two of the magnets to create a magnet slightly larger than the one on the edge of the face.
5. Use the north pole of the larger magnet to push the north pole of the smaller magnet onto the face. Explain that the Galvanic Spa Instrument uses the polarities or charges of particles with the charge of the head of the instrument to transport ingredients into the skin.
6. Now, if you place the north pole of the larger magnet toward the south pole of the smaller magnet, it will pull it out of the face. Explain that this represents how the instrument can pull toxins and impurities out of the skin with opposite polarities or charges.

FYI

Not every person needs the same level of current. By having different current levels (indicated by one, two, or three beeps), the Galvanic Spa Instrument allows each user to have the most comfortable and effective experience possible. To correctly adjust the current level to every individual, the instrument sends a controlled charge through the skin and measures how long it takes for the charge to return to the instrument. The amount of time it takes differs according to the moisture content, electrolyte balance, and overall health of the skin. The instrument then selects the level of current that will be most comfortable for your skin. This simple process replaces the years of training estheticians go through to effectively perform galvanic treatments in a spa setting.

Complementary Products

- **Galvanic Spa Body Shaping Gel**—works exclusively with the Galvanic Spa Instrument and Body Conductor to target pockets of fat and toxins for a firmer, more toned appearance.
- **Nutriol® Hair Fitness Treatment and Shampoo**—developed with patented technology, these advanced formulas remineralize the scalp and hydrate the hair with Tricalgoyl,® a clinically proven seaweed derivative. Protect the future of your hair with Nutriol.
- **Galvanic Spa Pre-Treatment and Treatment Facial Gels**—when your skin looks tired or stressed, a spa facial provides the revitalization you need.
- **Tru Face® Line Corrector**—in the ongoing search for key ingredients that may be paired with our exclusive Galvanic Spa technology, the Nu Skin Enterprises Research and Development team has discovered a simple, but powerful connection with Tru Face Line Corrector. When used with the positive current (two-fountain) setting on the Galvanic Spa instrument, the positively charged peptide in Tru Face Line Corrector more effectively permeates the skin, helping to deliver enhanced results.
- **Dermatic Effects® Body Contouring Lotion**—with hibiscus and malvaceae extracts that help inhibit fat production and stimulate fat breakdown, this contouring lotion aids in preventing the appearance of cellulite and aging skin. An eight-week clinical study showed a 55 percent reduction in the appearance of cellulite, an 84 percent

increase in moisture content, a 37 percent improvement in skin firmness, and a 69 percent enhancement in skin smoothness.

Frequently Asked Questions

How does the Galvanic Spa System differ from the new Galvanic Spa System II?

The new Galvanic Spa System offers more flexibility than the old Galvanic Spa System. With the Galvanic Spa System II or II EX, you can choose not only the polarity (charge) of a treatment, but also the length of time and the attachment. This allows you to easily customize your home spa experience.

How long will the batteries last?

In-house testing showed that batteries in the Galvanic Spa instrument lasted several months with normal usage.

Why does the Galvanic Spa instrument beep?

When you begin a treatment, you will hear one, two, or three beeps as the instrument automatically sets the current level for your skin type. You will also hear beeps at 20 second intervals during your treatment. These beeps were added to let the user know the treatment is in progress because the instrument was designed to be silent when in use. In addition, you will hear two beeps and a check mark will appear in the instrument's display when your treatment is complete.

The new instruments feel lighter than the last. Why is that?

Because Galvanic Spa II EX treatments can last up to five minutes, the instrument was designed to be lighter to reduce muscle fatigue.

How do I change the conductor attachments?

To release a conductor: With your free hand on the conductor being released, hold the instrument with the chrome activation panel facing you and press on the release button firmly with your thumb until the conductor is released from the instrument.

To attach a conductor: Hold the instrument with the display side up and hold the conductor with the positioning indicator (small etch) facing you. Press firmly to snap the conductor into place. If the conductor fails to attach, make certain the conductor is positioned correctly. Forcing the conductor could damage both the instrument and the conductor. The instrument will not work unless the conductor is attached correctly.

To reduce the risk to your health: Check with your physician before using the Galvanic Spa instrument or receiving a Galvanic Spa treatment if you are pregnant, have a pacemaker, are epileptic, wear metal braces on your teeth, or have a metal implant. The Galvanic Spa instrument should not be used over an area with many broken capillaries or over an inflamed or open wound.

* Use product only for its intended use as described in the Galvanic Spa System Owner's Manual.