



光感素肌 一觸即現  
TOUCH OF GLOW





作為 2017 至 2023 年全球排名第一美容儀器系列品牌\*  
，Nu Skin 永不停步，結合領先科技及優質成分，研  
發優質、高效及創新的產品，而 2020 年的重量級創  
新產品便是 ageLOC Boost™。

As the World's #1 brand for Beauty Device Systems from 2017  
to 2023\*, Nu Skin never stop to get a step further. Combining  
advanced technology and premium ingredients, Nu Skin  
continues to develop high quality and innovative products,  
and the spotlight of 2020 goes to our brand new ageLOC  
Boost™.

\*資料來源：歐睿國際有限公司；零售價值 RSP 條款、所有銷售管道；2017 年至 2022 年。美容系列指的是專屬搭配或建議與同品牌外用消費品一起使用之居家護膚美容儀器。宣稱驗證根據 2023 年 4 月至 6 月間所做的客製化研究與方法調查。居家護膚美容儀器銷售額包含Passport資料庫所定義之電動潔面儀的銷售額。此類別不包含護髮/除毛用具、剃鬚刀、和口腔保健用具。

\*Source: Euromonitor International Limited; Retail Value RSP terms; all channels; 2017 to 2022. Beauty Systems are at-home Skin Care Beauty Devices that are exclusively paired or recommended to be used with a topical consumable of the same brand. Claim verification based on Euromonitor custom research and methodology conducted April -June of 2023. Sales of at-home skin care beauty devices includes sales of electric facial cleansers as defined in Passport database. This category does not include hair care/removal appliances, body shavers, and oral care appliances.



## ageLOC boost™

每日亮澤

**Glow Everyday**

2 分鐘煥發飽滿光澤，塑造年輕彈力肌  
Optimizes and awakens skin's natural beauty  
for a youthful, glowing complexion in just 2  
minutes



## ageLOC LumiSpa® iO

早晚煥膚

**Renew Morning & Evening**

淨肌煥膚，一機雙效  
配搭睛亮修護組合讓雙眸綻放睛采  
Dual action for skin cleansing and renewal  
Use with Accent for an invigorating treatment  
to skin around the eyes



## 極致 Home Spa 體驗 Ultimate Home Spa Experience

Nu Skin 居家美容儀  
全面覆蓋您的護膚需求

Nu Skin device platform  
fully covered your skincare needs



## ageLOC® 修身美顏 Spa 及美體 Spa\* ageLOC® Galvanic Spa® & Body Spa®\*

每周護養

**Nourish Weekly**

促進主要成分輸送至肌膚  
一機多效，集多項修護抗衰老功能於一身  
Enhances the delivery of key ingredients to the skin  
Offers a wide range of anti-aging treatments in one  
single device



## ageLOC Me® 新智我® ageLOC Me®

日夜修護

**Treat Day & Night**

打造客製化抗衰老護膚方案，提供每日  
專屬修護  
Customized anti-aging skincare solution for  
everyday treatment



## ageLOC boost™

符合人體工學，舒適護膚體驗

Ergonomic design for comfortable treatment experience

輕巧便攜，隨時隨地煥活肌膚

Handy for a touch of glow anytime, anywhere

無按鍵設計，感應式充電，型格易用

Smart, fancy buttonless design and magnetic charger



隨時隨地 · 光采自生  
GLOW FROM WITHIN  
ANYTIME · ANYWHERE

Nu Skin 全新 ageLOC boost™ 以獨有「變向脈衝微電流」技術煥活肌膚。正負極均衡的脈衝每 2.4 秒由正脈衝切換成負脈衝，每 2 分鐘約 1000 次的脈衝讓肌膚告別惰性，帶來新鮮活膚體驗。當與光感瞬活精華同時使用，更能煥活年輕光采肌膚。

Nu Skin's newest innovation - ageLOC boost™ - uses unique "Variable Pulse Microcurrent" technology to activate the skin. The pulses are charge-balance and switch from positive to negative every 2.4 seconds. Approximately 1000 times microcurrent offers a fresh beauty boost to the skin in just 2 minutes! When used with Activating Treatment, the microcurrent optimizes and awakens skin's natural beauty for a youthful, glowing complexion.





ageLOC boost™

## 光感瞬活精華 Activating Treatment

配合 ageLOC boost™ 使用令肌膚更亮澤、飽滿及富彈性

---

Use with ageLOC boost™ for a brighter, plumper and bouncier skin

光感亮澤配方 — 蘊含於可控環境農業 (CEA) 種植的海百合，有效亮澤、均勻膚色

---

Intensive Whitening Blend – With *Pancratium Maritimum* grown through Controlled Environment Agriculture (CEA) for a bright, even skin tone

瞬活彈潤配方 — 讓肌膚水潤柔軟，塑造年輕飽滿彈力

---

Plumping Blend – Hydrates and soften for a fuller, plumper appearance

ageLOC® 專有配方 — 針對初老徵狀的根源，煥發年輕光采肌膚

---

ageLOC® Proprietary Blend – Target early signs of aging for youthful looking skin




# 實驗證實 12 星期顯著美肌

## CLINICALLY PROVEN BEAUTY BENEFITS IN 12 WEEKS


實驗證明，測試者在連續使用 ageLOC® Boost™ 12 星期後，肌膚於平滑、柔軟、亮澤、緊緻及膚色均勻的改善均持續提升#。

Clinical proven improvements on skin smoothness, softness, luminosity, firmness and skin tone evenness are continued to grow when used for 12 weeks#.



而測試者在使用 12 星期後的自我評分中，更感到肌膚彈性、飽度、毛孔顯著度及水潤度得到改善#。

Users even reported the improvements in skin bounce resilience, plumpness, pore visibility and hydration after used for 12 weeks#.



#Nu Skin (2022). A clinical study to evaluate the skin benefits of using ageLOC Boost™ with ageLOC Boost™ Activating Treatment on 25-40-year-old Asian females. Unpublished data.






# 使用方法 HOW TO USE

1/



用手指於輕觸式面板由下而上掃過啟動儀器。  
Turn on device by swiping finger up full length of touch surface.

2/



按壓 3 泵光感瞬活精華在 ageLOC boost™ 導頭上，均勻塗抹於半邊面上。  
Dispense 3 pumps of Activating Treatment onto the treatment surface. Spread over one half of the face.

3/



用手指由下而上掃過輕觸式面板，然後於面上以向上及向外的方式進行 1 分鐘護理。  
Begin treatment by swiping finger up full length of touch surface. Glide treatment surface over face in an upward / outward motion for 1 minute treatment.

4/



1 分鐘護理後儀器會暫停，請於另一邊面重覆步驟 2 及 3。按摩面上剩餘的精華直至吸收，然後塗上後續精華液及乳霜。  
When the device pause after 1 minute, repeat STEP 2 & STEP 3 on other half of the face. Massage excess Activating Treatment into skin, then follow with daily serum and moisturizer.

## Nu Skin 全方位護膚方案 360° SKINCARE SOLUTION

