

# **ageLOC®** **GALVANIC BODY SPA™**

DELIVERS UP TO 10 TIMES MORE AGELOC INGREDIENTS TO YOUR SKIN



## AGELOC POSITIONING

Nu Skin's exclusive ageLOC science targets arSuperMarkers—the ultimate sources of aging.

## PRODUCT OVERVIEW

*Your body optimized.* Nu Skin has captured the future of body care with the ageLOC® Galvanic Body Spa™. Ergonomically designed for ease of use and handling during body treatments, the new ageLOC Galvanic Body Spa with pulsating technology features a proprietary ageLOC body conductive surface designed to deliver 10 times more ageLOC to your skin when used with ageLOC Body Shaping Gel. The ageLOC Galvanic Body Spa targets the ultimate sources of aging while diminishing the appearance of fat and cellulite and smoothing skin to provide a slimmer, more toned appearance. It also features a self-adjusting current, lighted display, and audio indications, making it an easy to use in-home anti-aging body system.

Newly developed, patent-pending galvanic pulsating technology is optimized to maximize anti-aging effects on the arms, abdomen, buttocks, and thighs. The skin on the body is thicker and less sensitive, making it benefit from and respond to a pulsating galvanic current. It also helps stimulate, purify, and refresh your skin to reduce the visible signs of aging—for a slimmer, smoother, firmer, youthful looking you today and in the future.

## TARGET AUDIENCE

### PRIMARY MARKETS

- Women and men looking for a firmed, toned, younger looking body.
- Current ageLOC Galvanic Body Spa and Galvanic Spa® customers.

### SECONDARY MARKETS

- Image conscious women, ages 25–65, who may occasionally visit day spas.
- Women of any age with cellulite and other signs of aging on the body.

## PRODUCT BENEFITS

- Aids in delivering key ingredients to the skin.
- Optimized pulsating galvanic current aids in stimulating skin, helps purify and refresh skin, and helps reduce the visible signs of aging.
- Laboratory research demonstrates that treatments with a Galvanic current can enhance the delivery of anti-aging ingredients for up to 24 hours.

## ENHANCED FEATURES

- **Automatically adjusts to your skin**—delivers the correct amount of current for a more effective treatment.
- **Enhanced ergonomics**—the new shape of the ageLOC® Galvanic Body Spa™ is designed to better follow the contours of the body.
- **ageLOC conductive surface**—features a larger proprietary ageLOC conductive surface, which is able to maintain 10 times more product between the skin and the instrument surface, delivering more key ageLOC ingredients to the skin.

## USAGE

For optimal results, use your ageLOC Galvanic Body Spa with ageLOC Body Shaping Gel once a day three days a week on thighs, arms, buttocks, and abdomen.

Read through the owner's manual completely before starting your galvanic treatment.

## COMPLEMENTARY PRODUCTS

- **ageLOC® Body Shaping Gel**—formulated to work exclusively with Nu Skin® Galvanic Spa® instruments to bring ageLOC anti-aging benefits to the body, smooth the appearance of fat and cellulite, and help refresh and purify the skin while reducing the visible signs of aging.
- **ageLOC® Dermatic Effects Body Contouring Lotion**—scientifically formulated to bring ageLOC anti-aging benefits to the body. This daily moisturizing and contouring lotion contains ingredients that help inhibit fat production, stimulate fat breakdown, and prevent the appearance of fat and cellulite while smoothing the skin.

## FREQUENTLY ASKED QUESTIONS

**What is unique about the ageLOC Galvanic Body Spa?**

The ageLOC Galvanic Body Spa was designed for use in combination with ageLOC Body Shaping Gel on body areas—not on the face and neck. This instrument was designed to target a larger surface area with a patent-pending ageLOC conductive surface that captures 10 times more key ingredients between the conductive surface and the skin (when compared to the body conductor attachment). Secondly, this design is more ergonomic for easier use and handling during body treatments. Thirdly, the ageLOC Galvanic Body Spa and ageLOC Body Shaping Gel target end benefits/results (i.e., appearance of cellulite, slimming, toning, etc.) that are different from those targeted on the face. Lastly, the skin on the body is thicker and less sensitive, making it more apt to benefit from and respond to a pulsating current.

**What is the difference between a pulsating galvanic current and a constant galvanic current?**

A constant galvanic current is one that has a steady flow. It does not waver or vary in amplitude. A pulsating galvanic current is one that varies between two different, pre-established levels. Please note these are very different than alternating currents (AC). Within alternating currents, the polarity fluctuates between positive and negative poles. A pulsating galvanic current is better suited for the body.

**What is the benefit of using a pulsating galvanic current on the body?**

Direct (galvanic) currents, whether pulsating or not, can assist in the delivery of like charged key ingredients to the skin by improving osmotic flow. This is increasingly important for enhancing delivery of key ingredients to areas with thicker skin. Additional research has shown that the pulsating current enhances fluid vibration. Pulsating galvanic currents aid in stimulating skin and increasing cellular communication to help purify and refresh skin and help reduce the visible signs of skin aging. Pulsating galvanic currents are optimally designed for use on the body where the skin can benefit from enhanced fluid vibration.

**What is the benefit of having a grooved body conductive surface?**

The patent-pending grooved ageLOC conductive surface captures 10 times more key ingredients in the critical area between the conductive surface and the skin (when compared to the body conductor attachment).

**How much pressure should I apply?**

Nu Skin's Galvanic Spa instruments should be very gently used with light pressure. There is no need to apply pressure to get the maximum benefits. The galvanic currents do the work for you.

**Should the ageLOC Galvanic Body Spa be used on the face?**

No. The ageLOC Galvanic Body Spa was not designed for use on the face and neck. This instrument was designed specifically for use on the body. It easily treats a larger surface area, easily follows the contours of the body, and, when combined with ageLOC Body Shaping Gel, brings about the targeted end results (i.e., diminishing appearance of cellulite, smoothing skin, slimming, toning, etc.). It also provides a pulsating current, which is better suited for structural differences and targeted attributes in the area being treated. The ageLOC® Galvanic Body Spa™ should only be used as instructed.

**Are the ageLOC Galvanic Body Spa, ageLOC Body Shaping Gel, and ageLOC Dermatic Effects weight-loss products?**

No, the ageLOC Galvanic Body Spa, ageLOC Body Shaping Gel, and ageLOC Dermatic Effects are not

# ageLOC® GALVANIC BODY SPA™

weight-loss products. ageLOC Body products were developed to minimize and smooth the appearance of fat and cellulite, improve the appearance of skin's firmness, and promote the youthful appearance of the skin.

**TO REDUCE THE RISK TO YOUR HEALTH:** The ageLOC Galvanic Body Spa instrument is only for use on healthy skin. Do not use if you have open sores or wounds, skin inflammation, allergies to metal, or problematic or overly sensitive skin. Check with your physician before using if you are pregnant, have a pacemaker or similar device, are epileptic, have any metal implants, or in case of illness.