

ageLOC™ EDITION GALVANIC SPA® SYSTEM II

BOOSTS VISIBLE BENEFITS OF AGELOC™
FUTURE SERUM BY AN AVERAGE OF 80%*

ageLOC™ POSITIONING

Nu Skin's exclusive ageLOC science targets arSuperMarkers—the ultimate sources of aging.

PRODUCT OVERVIEW

Nu Skin has unlocked the future of skin care with the new, revolutionary ageLOC Edition Galvanic Spa System II and ageLOC Face Conductor. This powerful combination is designed to deliver up to **five times more ageLOC ingredients** to the skin when used with Galvanic Spa Facial Gels with ageLOC⁺—targeting the ultimate sources of aging.

Now you can enjoy the benefits of a spa visit as a part of your regular anti-aging skin care regimen. The small, smart, and powerful ageLOC Edition Galvanic Spa System II features patented self-adjusting currents and interchangeable conductors that work synergistically with Nu Skin's specially formulated treatment products to facilitate the transport of key ingredients to the skin. The newly engineered ageLOC Galvanic Spa also includes a larger, brighter display, intuitive usability, and enhanced ergonomics—making it the easiest to use, most cost effective in-home spa system available.

Now you can unlock the secret to a younger looking you with the ageLOC Edition Galvanic Spa System II. Your radiant, youthful looking face and body will leave people wondering about your age—and your secret.

BENEFITS

- Laboratory research demonstrates that treatments with Galvanic current can enhance the delivery of anti-aging ingredients for up to 24 hours.

* When using ageLOC Future Serum daily, the central product of ageLOC Transformation, and the ageLOC Galvanic Spa three times a week, as determined by clinical assessment after seven days of use.

† When compared to the Galvanic Spa II instrument.



- The Galvanic Spa instrument features proprietary, patented technology with four interchangeable conductors:
 - **ageLOC Face Conductor**—the larger contoured and grooved surface area keeps product in contact with the skin and provides the skin with up to five times more ageLOC.
 - **Spot Treatment Conductor**—a flat, smooth, oval-shaped treatment node centralizes surface contact area, providing better contact with the skin, helping deliver more concentrated treatment to specific areas.
 - **Body Conductor**—three rounded nodes maximize body surface area contact for complete coverage and better massaging action.
 - **Scalp Conductor**—the comb-like conductor easily glides through the hair, while maintaining crucial contact between the conductor and your scalp.

ENHANCED FEATURES

- **Intuitive usability**—no more guesswork. The timing and polarity are automatically pre-programmed to adjust accordingly to your treatment selection.
- **Enhanced ergonomics**—the new shape of the ageLOC Face Conductor is designed to better follow the contours of the face, while the grooved surface keeps the facial gels in constant contact with your skin.
- **Automatically adjusts to your skin**—delivers the correct amount of current for a more effective treatment.
- **Lighted display**—new, easy-to-read LCD display includes a backlit icon to inform you that the instrument is in active mode.
- **Audio indications**—beeps every 10 seconds and has new “end treatment” and “low battery” notification.

ageLOC™ EDITION GALVANIC SPA® SYSTEM II

DEMONSTRATION

Before performing this magnet demonstration, explain that all galvanic instruments use electrical charges to facilitate the transport of key ingredients for maximum benefit. These instruments utilize an important characteristic of electrical charges—like charges repel each other and opposite charges attract each other.

(+) ◀▶ (+) (-) ◀▶ (-)

(+) ▶◀ (-) (-) ▶◀ (+)

If a key ingredient is negatively charged, the negative cycle of the Galvanic Spa Instrument will repel this ingredient—facilitating its transport. In the same way, if the ingredient is positively charged, the instrument will facilitate transport when the cycle is positive. If there are impurities within the skin that need to be pulled out, they will be extracted with an opposite charge.

To demonstrate how Galvanic Spa System currents work synergistically with specially formulated products to facilitate the transport of select key ingredients, use magnets on a magnetic board. Follow the steps below:

1. Obtain a magnetized board and several magnets (it's best if the magnets are wheel-shaped so they will roll across the board).
2. Draw or paste a picture of a face onto the magnet board.
3. Place one of the wheel-shaped magnets on the edge of the face.
4. Connect two of the magnets to create a magnet slightly larger than the one on the edge of the face.
5. Use the north pole of the larger magnet to push the north pole of the smaller magnet onto the face. Explain that the Galvanic Spa Instrument uses the polarities or charges of particles with the charge of the head of the instrument to transport ingredients into the skin.
6. Now, if you place the north pole of the larger magnet toward the south pole of the smaller magnet, it will pull it out of the face. Explain that this represents how the instrument can cleanse toxins and impurities from the skin with opposite polarities or charges.

FYI

Not every person needs the same level of current. By having different current levels, the Galvanic Spa Instrument allows each user to have the most comfortable and effective experience possible. To correctly adjust the current level to every individual, the instrument sends a controlled charge through the skin and measures how long it takes for the charge to

return to the instrument. The instrument then selects the level of current that will be most comfortable for your skin. This simple process replaces the years of training estheticians go through to effectively perform galvanic treatments in a spa setting.

GALVANIC SPA COMPATIBLE PRODUCTS

- **Galvanic Spa® Facial Gels with ageLOC™**—target the ultimate sources of aging, what Nu Skin has coined “arSuperMarkers”; for immediate and long-term benefits. Help stressed, tired skin immediately look younger and revitalized.
- **Tru Face® Line Corrector**—helps to decrease the appearance of fine lines and wrinkles. When used with the Galvanic Spa, Tru Face Line Corrector provides faster and better results than when used alone.
- **Galvanic Spa® Body Shaping Gel**—works exclusively with the Galvanic Spa Instrument and Body Conductor to target pockets of fat and toxins for a firmer, more toned appearance.
- **Nutriol® Hair Fitness Treatment**—developed with patented technology, this advanced formula works with the Scalp Conductor to help hair look and feel more abundant. Features a clinically proven seaweed derivative, Tricalgoxyl®, and is formulated to re-mineralize the scalp and hydrate the hair.

CONSUMER USE STUDY

- 100% of users in an in-house consumer use study indicated that the ageLOC Edition Galvanic Spa System II had “very good” to “excellent” overall readability of the display.
- 96% of users in an in-house consumer use study agreed that selecting treatment with the ageLOC Edition Galvanic Spa System II is simple.
- 91% of users in an in-house consumer use study graded the ergonomics of the ageLOC Edition Galvanic Spa System II as “good” to “excellent.”
- 91% of users in an in-house consumer use study graded the effectiveness of the new ageLOC Face Conductor as “good” to “excellent.”

TESTIMONIALS

- “[The] ageLOC [conductor] was a great feeling and it seems to fit with the natural curve of my face.”—Ali F.
- “I loved the display on the ageLOC unit; it was so much brighter.”—Lisa H.
- “I love the ageLOC unit!”—Kim C.

ageLOC™ EDITION GALVANIC SPA® SYSTEM II

- “The ageLOC unit’s new head shape was great!”
—Zachary L.
- “I liked the feel of the conductor of the ageLOC unit; it was easier to cover [the] whole side of my face; I also liked the clock display.” —Christy L.

COMPLEMENTARY PRODUCTS

- **ageLOC™ Elements**—three premier products target aging at its source with ageLOC while working to cleanse, purify, renew, moisturize, and reveal a radiant, younger looking complexion.
- **ageLOC™ Transformation**—infuse your skin with the power of ageLOC™ on a daily basis to target aging at its source. These four premier products work together to reveal younger looking skin in eight ways—for a more youthful, healthier looking you now and in the future.
- **ageLOC™ Future Serum**—delivers the maximum concentration of ageLOC to target the ultimate sources of aging. This luxurious serum is clinically proven to reveal younger looking skin in eight significant ways: by promoting youthful skin structure, smoothing texture, reducing the appearance of fine lines and wrinkles, increasing radiance, minimizing the appearance of pores, evening skin tone, reducing discoloration, and increasing hydration.
- **Tru Face® Essence Ultra**—restores skin definition and firmness. Features Ethocyn®, an anti-aging ingredient that has been clinically proven to improve the skin’s elastin content, the key component to firm skin.
- **Dermatic Effects® Body Contouring Lotion**—with hibiscus and malvaceae extracts that help inhibit fat production and stimulate fat breakdown, this contouring lotion aids in preventing the appearance of cellulite and aging skin.
- **Nutriol® Shampoo**—features Tricalgoxyl®, a clinically proven seaweed derivative that works with other essential ingredients to re-mineralize the scalp and help hair regain its natural health, strength, and shine.

FREQUENTLY ASKED QUESTIONS

What is the benefit of having a grooved Facial Conductor?

The grooves on the ageLOC Edition Galvanic Facial Conductor create an area that captures more ageLOC product between the conductor and the skin. With additional product in this critical location, the delivery of key ageLOC ingredients to the skin is enhanced. As the ageLOC Facial Conductor is moved around the skin, fresh product is captured and moved into the critical area directly between the skin and conductor. This helps maintain a more

consistent concentration of ageLOC actives on the skin.

How long will the batteries last?

In-house testing showed that batteries in the Galvanic Spa instrument lasted several months with normal usage.

Why does the Galvanic Spa instrument beep?

When you begin a treatment, you will hear one, two, or three beeps as the instrument automatically sets the current level for your skin type. You will also hear beeps at 10-second intervals during your treatment. In addition, you will hear one long beep followed by two shorter beeps when your treatment is complete.

How do I change the conductor attachments?

To release a conductor: With your free hand on the conductor being released, hold the instrument with the chrome activation panel facing you and press on the release button firmly with your thumb until the conductor is released from the instrument.

To attach a conductor: Hold the instrument with the display side up and hold the conductor with the positioning indicator (small etch) facing you. Press firmly to snap the conductor into place. If the conductor fails to attach, make certain the conductor is positioned correctly. Forcing the conductor could damage both the instrument and the conductor. The instrument will not work unless the conductor is attached correctly.

To reduce the risk to your health: Check with your physician before using the Galvanic Spa instrument or receiving a Galvanic Spa treatment if you are pregnant, have a pacemaker or similar device, are epileptic, wear metal braces on your teeth, or have a metal implant. The Galvanic Spa instrument should not be used over an area with many broken capillaries or over an inflamed or open wound. Use product only for its intended use as described in the Galvanic Spa System Owner’s Manual.