




SYSTEM OVERVIEW

Are you ready for a first-of-its-kind, dual-action skin care system that delivers dramatic skin renewal and deep-cleansing in a single treatment? Developed in conjunction with a board-certified dermatologist, ageLOC LumiSpa is inspired by gene expression studies of skin protein production.

With its soft silicone treatment head, ageLOC LumiSpa moves at the precise frequency to promote skin renewal, revealing smoother, softer skin after just one use.

LumiSpa gently—yet deeply—cleanses away dirt, oil, makeup, and toxins through its patent-pending counter-rotating, pore-tightening action. It leaves skin smoother and purer, while reducing the appearance of pores. And the stimulating, massaging LumiSpa experience will leave you looking forward to your next treatment.

Experience brighter, softer, smoother skin and an energizing facial massage. Simply spend two minutes twice a day to enjoy healthier, youthful looking skin in as little as two weeks.

TARGET AUDIENCE

ageLOC LumiSpa is suitable for any skin-conscious adult seeking ways to minimize the visible signs of aging or to maintain and promote healthy, youthful looking skin.

KEY FEATURES

- **Device**—ageLOC LumiSpa is a waterproof and rechargeable handheld personal skin care device. Used together with the specially formulated LumiSpa treatment cleansers, it simultaneously treats and cleanses the face.

- **Treatment Heads**—ageLOC LumiSpa features two treatment heads: **NORMAL** and **FIRM**. Their gentle, non-abrasive silicone surfaces contains antimicrobial silver and is more hygienic than brushes commonly used with other devices. Their sophisticated design helps to maximize product efficacy while making skin more receptive to subsequent products applied to the skin.
- **Treatment Cleansers**—ageLOC LumiSpa features four treatment cleanser options designed for different skin types: dry skin, normal/combo skin, oily skin, and sensitive skin. These powerful products are formulated using Nu Skin's proprietary ageLOC ingredient blend. Optimally designed for LumiSpa, they are proven to maximize the skin perfecting and cleansing benefits of the device—revealing healthy, bright, youthful looking skin.

BENEFITS

- Provides instant benefits so skin feels softer and smoother after just one use.
- Delivers seven skin benefits: skin renewal, smoothness, radiance, a more refreshed look, purified skin, minimized appearance of pores, and increased appearance of volume and density.
- LumiSpa benefits intensify over time when consistently used for two minutes twice daily.
- When used in the morning, LumiSpa energizes the skin for a fresher, smoother, more rejuvenated complexion.
- When used at night, LumiSpa deeply cleanses the skin, helping massage and de-stress the skin.
- Removes dead skin cells, promoting dynamic skin renewal and improved skin brightness and smoothness.
- Soft silicone treatment heads are more hygienic, less abrasive, and easier to clean than bristle heads commonly found on the market.

- Delivers brighter, healthier, more youthful looking skin in as little as two weeks.
- Gentle enough for twice-daily use.
- Unique sweeping action lifts away loose cells while reducing appearance of pores.
- Goes beyond cleansing to deliver a skin-perfecting treatment through its patented counter-rotating design.
- Reduces appearance of visible pores.
- Doesn't strip skin, but retains skin's natural moisture balance as it cleanses.
- Gently cleanses away impurities, makeup, dirt, and toxins without irritating skin.
- Formulated with Nu Skin's proprietary ageLOC ingredient blend to target the sources of aging and preserve the look of youth.

CLINICALLY PROVEN BENEFITS

A preliminary 12-week study using ageLOC LumiSpa revealed that individuals who used the product experienced improvements in skin:

- Smoothness
- Softness
- Texture
- Clarity
- Radiance

Participants noticed several of these improvements after a single use.

USAGE

For optimal results, use ageLOC LumiSpa for two minutes twice daily, as your cleansing step. Morning and night, treat and cleanse with ageLOC LumiSpa treatment cleanser. Then follow with your preferred toner, serum, and moisturizer.

AGELOC® LUMISPA™ TREATMENT CLEANSERS

ageLOC LumiSpa features four treatment cleanser options designed for different skin types: dry, normal/combo, oily, and sensitive skin. These powerful products contain Nu Skin's proprietary ageLOC ingredient blend to target the sources of aging and are formulated to maintain healthy, bright, youthful looking skin.

ageLOC treatment cleansers are specifically formulated to enable precise cushioning, cleansing, and interaction with the skin, optimizing the effectiveness of LumiSpa's patent-pending treatment head action.

AGELOC® LUMISPA™ TREATMENT CLEANSERS KEY INGREDIENTS

- **ageLOC LumiSpa Treatment Cleanser—Oily** contains Japanese Knotweed (Polygonum Cuspidatum Root) Extract, Myristyl Alcohol, and PCA to help reduce oil production and soothe the skin, as well as Pomegranate (Punica Granatum) Juice Extract and Carnosine to provide antioxidant protection.
- **ageLOC LumiSpa Treatment Cleanser—Normal/Combo** contains Milk Thistle (Silybum Marianum Seed) Oil, Porphyra Umbilicalis Extract, and Nannochloropsis Oculata Extract to help maintain skin's natural moisture barrier, as well as Rhododendron (Ferrugineum) Extract to provide antioxidant protection.
- **ageLOC LumiSpa Treatment Cleanser—Dry** contains Sodium Lactate, Sodium PCA, Glycine, Fructose, Niacinamide, and Inositol to help protect skin's natural moisture, Squalane to help maintain skin's natural moisture barrier, and Salicyloyl Phytosphingosine to soothe the skin.
- **ageLOC LumiSpa Treatment Cleanser—Sensitive** is our most gentle formula. It is a pink gel cream formulated with a mild surfactant system, which results in a low foam. It contains Avena Sativa (Oat) Kernel Extract, Pinus Tabulaeformis Bark Extract (which contributes to the pink color), and Bisabolol to help calm and soothe the delicate skin. Allantoin also helps to condition the skin.

COMPLEMENTARY PRODUCTS

- **ageLOC Galvanic Spa/Nu Skin Facial Spa**—The ageLOC Galvanic Spa utilizes direct/galvanic currents to rejuvenate your complexion and increase skin vibrancy. The Nu Skin Facial Spa uses microcurrents to stimulate and tone the skin. Use Galvanic Spa or Facial Spa two to three times weekly along with daily ageLOC LumiSpa treatments to dramatically increase visible results.
- **ageLOC Me®**—Creates an entirely new approach to customized anti-aging skin care. Brilliant yet simple to use, ageLOC Me delivers five powerful anti-aging products, custom selected from almost 2,000 possible regimen combinations. Delivering one streamlined, simplified regimen to you.
- **Tru Face® Line Corrector**—Targets signs of aging with peptides that help soften the appearance of lines around your mouth, eyes, and forehead.
- **ageLOC® Tru Face® Essence Ultra**—Formulated with the power of Ethocyn®* and Nu Skin's revolutionary anti-aging science to target the sources of aging that lead to the loss of firmness. Leave undefined skin in the past and

enjoy firmer, younger looking skin today and in the future with ageLOC Tru Face Essence Ultra.

FREQUENTLY ASKED QUESTIONS

How does ageLOC LumiSpa work?

With its soft silicone treatment head, ageLOC LumiSpa's proprietary Micropulse Oscillation technology enables the device to move at the precise frequency to promote skin renewal, revealing smoother, softer skin after just one use. At the same time, LumiSpa gently—yet deeply—cleanses away dirt, oil, makeup, and toxins through its patented counter-rotating, pore-tightening action, leaving skin smoother and purer while reducing the appearance of pores.

What are the main differentiators of ageLOC LumiSpa?

There are three key differentiators of ageLOC LumiSpa:

1. Unlike cleansing devices currently on the market, which feature abrasive brush or bristle heads, ageLOC LumiSpa features a gentle, more hygienic silicone treatment head that won't damage the skin or lift skin cells.
2. ageLOC LumiSpa features patented Micropulse Oscillation technology to stimulate the skin to more effectively treat and cleanse the skin and deliver superior skin renewal benefits.
3. ageLOC LumiSpa is a dual-action, timesaving device that delivers a superior skin-perfecting treatment and deep cleansing benefits simultaneously.

What is the rationale behind the name?

"Lumi" stands for luminous and emphasizes the initial benefits that using the device and products will have on your skin. "Spa" reinforces the fact that the device, treatment heads, and treatment cleansers will deliver spa-quality results every time you use them.

How do I know which ageLOC LumiSpa treatment head is right for me?

Select your ageLOC LumiSpa treatment head based on your skin type and treatment needs. Both treatment heads remove impurities and purify skin, delivering a gentle experience ideal for sensitive, dry, oily, and normal to combination skin. The Firm treatment head offers a more aggressive cleansing option.

How long will my treatment head last?

We recommend replacing your treatment head every three months to ensure your skin will receive optimal benefits on a consistent basis.

How does ageLOC LumiSpa work with the treatment cleansers?

The treatment cleansers were developed specifically for use with ageLOC LumiSpa. They're designed to work synergistically with the counter-rotating motion of the heads to optimize distribution of anti-aging ingredients onto the skin while cleansing and exfoliating to promote healthy, youthful looking skin.

How long should the treatment cleanser last?

One treatment cleanser should last 30 days.

How many studies has Nu Skin conducted for ageLOC LumiSpa?

Nu Skin conducted 11 studies during ageLOC LumiSpa development and will be conducting additional clinical studies to substantiate benefits in the near future.

Are there any patents on ageLOC LumiSpa?

Yes, there are multiple patents issued and pending for ageLOC LumiSpa and the proprietary Micropulse Oscillation technology.

Do I need to use the ageLOC LumiSpa treatment cleansers with ageLOC LumiSpa, or can I use other treatment formulas or cleansers with it?

ageLOC LumiSpa treatment cleansers are specifically formulated for use with the ageLOC LumiSpa device. They enable the precise cushioning, cleansing, and interaction with the skin, optimizing the effectiveness of LumiSpa's patented treatment head action. They are the only formulations tested and approved for use with ageLOC LumiSpa to deliver the full system benefits and promote healthy, youthful looking skin.

How is ageLOC LumiSpa different from ageLOC Galvanic Spa® or Nu Skin® Facial Spa?

ageLOC LumiSpa provides twice daily treatment and cleansing and should be used as part of your regular skin care routine, whereas ageLOC Galvanic Spa and Nu Skin Facial Spa are specialized treatments that can be used three times per week. ageLOC LumiSpa and ageLOC Galvanic Spa/Nu Skin Facial Spa deliver complementary anti-aging benefits: ageLOC LumiSpa treats and cleanses skin for a smoother, softer, more radiant appearance while ageLOC Galvanic Spa conditions, revives, and hydrates the skin, and Nu Skin Facial Spa uses microcurrents to tone and contour the skin.

Does ageLOC LumiSpa treat acne?

ageLOC LumiSpa was not designed to treat acne. The system was developed to promote healthy, youthful looking skin and to deliver general skin renewal and purification benefits. Those treating acne should continue with their regular skin care regimen and consult a dermatologist before using ageLOC LumiSpa.

Can I use ageLOC LumiSpa on my body?

ageLOC LumiSpa was designed to target facial skin and deliver anti-aging and deep cleansing benefits in this area. We do not recommend it be used elsewhere on the body.

Can I use ageLOC LumiSpa in the shower?

Yes, ageLOC LumiSpa is waterproof, and you can use it in the shower or other wet environments.

How long does ageLOC LumiSpa last after a full charge?

ageLOC LumiSpa was developed to last at least a week between charges with regular use.

4

INGREDIENTS

AGELOC LUMISPA TREATMENT CLEANSER—NORMAL/COMBO

Water (Aqua), Coco-Betaine, Sodium Cocoyl Glycinate, Glycerin, Acrylates Copolymer, Diisopropyl Sebacate, Glycol Distearate, PEG-8, Rhododendron Ferrugineum Extract, Nannochloropsis Oculata Extract, Schizandra Chinensis Fruit Extract, Porphyra Umbilicalis Extract, Silybum Marianum Seed Oil, Narcissus Tazetta Bulb Extract, Oryza Sativa (Rice) Extract, Acrylates/C10-30 Alkyl Acrylate Crosspolymer, Polyquaternium-39, Phenoxyethanol, Chlorphenesin, Sodium Benzoate, Phytic Acid, Sodium Hydroxide, Sodium Phytate, Fragrance (Parfum), Synthetic Fluorphlogopite, Titanium Dioxide (CI 77891), Alumina, Tin Oxide

AGELOC LUMISPA TREATMENT CLEANSER—OILY

Water (Aqua), Cocamidopropyl Betaine, Acrylates Copolymer, Ammonium Laureth Sulfate, Glycerin, PEG-80 Sorbitan Laurate, Di-PPG-2 Myreth-10 Adipate, Coco-Glucoside, Punica Granatum Juice Extract, Narcissus Tazetta Bulb Extract, Polygonum Cuspidatum Root Extract, Schizandra Chinensis Fruit Extract, Oryza Sativa (Rice) Extract, Glyceryl Oleate, Polyquaternium-10, Myristyl Alcohol, PCA, Carnosine, 1,2-Hexanediol, Fragrance (Parfum), Citric Acid, Aminomethyl Propanol, Dehydroacetic Acid, Chlorphenesin, Phenoxyethanol,

Sodium Benzoate, Propanediol, Caprylyl Glycol, Disodium EDTA, Mica (CI 77019), Titanium Dioxide (CI 77891)

AGELOC LUMISPA TREATMENT CLEANSER—DRY

Water (Aqua), Dicaprylyl Carbonate, Diisostearyl Malate, Tribehenin PEG-20 Esters, Glycerin, Glycereth-26, Squalane, Cetearyl Alcohol, Dimethicone, Caprylic/Capric/Myristic/Stearic Triglyceride, Sodium Astrocaryum Murumuru, Oryza Sativa (Rice) Extract, Schizandra Chinensis Fruit Extract, Xanthophyll, Narcissus Tazetta Bulb Extract, Acrylates/C10-30 Alkyl Acrylate Crosspolymer, Niacinamide, Urea, Methylsilanol Hydroxyproline Aspartate, PEG-40 Hydrogenated Castor Oil, Sodium Lactate, Sodium PCA, Vegetable Oil, Xanthan Gum, Batyl Alcohol, Butylene Glycol, Pentylene Glycol, Phytosteryl/Octyldodecyl Lauroyl Glutamate, Maltodextrin, Brassicamidopropyl Dimethylamine, Hydroxyphenyl Propamidobenzoic Acid, Salicyloyl Phytosphingosine, Fructose, Glycine, Inositol, Lactic Acid, Lecithin, Phenoxyethanol, Hydroxyacetophenone, Disodium EDTA, Sodium Hydroxide, Fragrance (Parfum)

AGELOC LUMISPA TREATMENT CLEANSER—SENSITIVE

Water (Aqua), Cetearyl Alcohol, Glycerin, PEG-8 Dimethicone, Cocamidopropyl Betaine, Caprylyl/Capryl Glucoside, Capryl/Capramidopropyl Betaine, Carbomer, Narcissus Tazetta Bulb Extract, Pinus Tabulaeformis Bark Extract, Schizandra Chinensis Fruit Extract, Avena Sativa (Oat) Kernel Extract, Oryza Sativa (Rice) Extract, Allantoin, Zingiber Officinale (Ginger) Root Extract, Bisabolol, Carnosine, Xanthan Gum, Hydroxyacetophenone, Phenoxyethanol, Disodium EDTA, Sodium Hydroxide, Sodium Chloride, Citric Acid