

ageLOC® LumiSpa®



MEET AGELOC LUMISPA

A first-of-its-kind, dual-action skin care system that helps to deliver dramatic skin renewal and deep cleansing in 1 Step for 2 Minutes with 7 Benefits.

When technology meets skin care, ageLOC LumiSpa is an innovation centered on a multi-award winning beauty device that works with soft silicone treatment head to move at the precise frequency needed to promote skin renewal—revealing softer, smoother skin after just one use.

ageLOC LumiSpa gently—yet deeply—cleanses away dirt, oil, makeup, and impurities through its proprietary oscillating, pore-purifying action. It leaves skin smoother and purer, while minimizing the appearance of pores.

Indulge in a rejuvenating facial beauty experience and enjoy softer, smoother skin that glows! Simply spend 2 minutes twice a day for a healthier and more youthful looking skin in as little as two weeks.


WHY YOU'LL LOVE IT

- In 1 Step, 2 Minutes, it helps to deliver 7 Skin Benefits: softness, smoothness, radiance, clarity, purified skin, minimized appearance of pores, and improved appearance of skin volume and density.
- Provides instant benefits so skin feels softer and smoother after just one use.
- Benefits increases over time when consistently used for 2 minutes twice daily.
- Removes dead skin, helps in promoting dynamic skin renewal and improved skin brightness and smoothness while retaining skin's natural moisture balance as it cleanses.

- Gentle enough to use twice daily. In the morning to refresh skin for a glowing, smoother, more rejuvenated complexion. At night to deeply cleanse the skin, helping to purify skin in a relaxing facial experience.

WHAT POWERS IT

- **Device**—ageLOC LumiSpa is a waterproof and rechargeable handheld personal skin care device which can be used in the shower or other wet environment.
- **Treatment Heads**—ageLOC LumiSpa features three treatment heads to personalize to your needs:



The soft, nonabrasive silicone surfaces are embedded with antimicrobial silver, more hygienic, less abrasive, and easier to clean than bristle heads commonly used with other skin care devices found in the market. The innovative designs help to maximize product efficacy while making skin more receptive to subsequent products.

- **Activating Cleansers**—ageLOC LumiSpa features five activating cleanser options designed for different skin types: normal to combination, oily, dry, sensitive and for blemish-prone. They are specifically formulated to enable the precise cushioning, cleansing, and interaction with the skin, optimizing the effectiveness of LumiSpa's proprietary treatment head action. These powerful products also contain Nu Skin's proprietary ageLOC ingredient blend for a more youthful look.

*Treat refers to the delivery of 7 skin benefits: softness, smoothness, radiance, clarity, purified skin, minimized appearance of pores, improved appearance of skin volume & density.
*ageLOC LumiSpa Activating Cleanser-For Blemish Prone is also known as ageLOC LumiSpa Treatment Cleanser-Acne in other markets outside of Southeast Asia.

ageLOC® LumiSpa®

Choose & Personalize To Your Skin



ageLOC LumiSpa Activating Cleanser—Normal/Combo contains Milk Thistle (Silybum Marianum Seed) Oil, Porphyra Umbilicalis Extract, and Nannochloropsis Oculata Extract to help maintain skin's natural moisture barrier, as well as Rhododendron (Ferrugineum) Extract to provide antioxidant protection.



ageLOC LumiSpa Activating Cleanser—Oily contains Japanese Knotweed (Polygonum Cuspidatum Root) Extract, Myristyl Alcohol, and PCA to help reduce oil production and soothe the skin, as well as Pomegranate (Punica Granatum) Juice Extract and Carnosine to provide antioxidant protection.



ageLOC LumiSpa Activating Cleanser—Dry contains Sodium Lactate, Sodium PCA, Glycine, Fructose, Niacinamide, and Inositol to help protect skin's natural moisture, Squalane to help maintain skin's natural moisture barrier, and Salicyloyl Phytosphingosine to soothe the skin.



ageLOC LumiSpa Activating Cleanser—Sensitive is our most gentle formula. It is a pink gel cream formulated with a mild surfactant system, which results in a low foam. It contains Avena Sativa (Oat) Kernel Extract, Pinus Tabulaeformis Bark Extract (which contributes to the pink color), and Bisabolol to help calm and soothe delicate skin. Allantoin also helps to condition the skin.



ageLOC LumiSpa Activating Cleanser—For Blemish Prone^{*} is a mild formula that contains 0.5% salicylic acid to help clear pores and reduce most breakouts. It is low-foaming to prevent additional skin irritation. It also contains Carnosine to provide antioxidant protection and Epilobium Angustifolium Flower/Leaf/Stem Extract that can help calm and support acne-prone skin. It is clinically proven to help mitigate acne breakouts and help the skin look smoother when used with ageLOC LumiSpa and Gentle treatment head.

DID YOU KNOW

ageLOC LumiSpa uses proprietary Micropulse Oscillation Technology, a pore-purifying action that draws out dirt, oil, makeup, and impurities and lifts them away without irritating skin while helping to support healthy skin proteins.

The unique sweeping action removes dead skin while reducing appearance of pores. It goes beyond cleansing to deliver a skin-perfecting experience through this proprietary oscillating design.

CLINICALLY PROVEN RESULTS

A 12-week study using ageLOC LumiSpa revealed that individuals who used the product experienced improvements in skin:

- Softness
- Smoothness
- Cleansing
- Brightness
- Radiance
- Clarity
- Luminosity
- Pore Appearance
- Firmness


Participants noticed several of these improvements after a single use.

HOW TO USE IT

For optimal results, use ageLOC LumiSpa for 2 minutes twice daily, as your cleansing step in the morning and at night with ageLOC LumiSpa Activating Cleanser. Then follow with your preferred toner, serum, and moisturizer.

Attach treatment head to the ageLOC LumiSpa device before you start:

1. Moisten face with water.
2. Apply ageLOC LumiSpa Activating Cleanser to all areas of the face.
3. Wet treatment head and press the power button to turn on device.
4. Slowly glide the treatment head back and forth over one of the areas of your face using broad strokes.



Tip: Do not use a scrubbing motion or press too hard.

5. Spend 30 seconds on each quadrant of your face. The device will pause quickly three times every 30 seconds to prompt you to move from one area of the face to the next.
6. When it's done, the device will stop and turn off automatically. Rinse your face with water to remove residual activating cleanser.

*Treat refers to the delivery of 7 skin benefits: softness, smoothness, radiance, clarity, purified skin, minimized appearance of pores, improved appearance of skin volume & density.
*ageLOC LumiSpa Activating Cleanser—For Blemish Prone is also known as ageLOC LumiSpa Treatment Cleanser—Acne in other markets outside of Southeast Asia.



TRY THESE WITH IT

- **ageLOC Galvanic Spa®**—uses direct/galvanic currents to rejuvenate your complexion and increase skin vibrancy.
- **ageLOC Me®**—creates an entirely new approach to customized anti-aging skin care. Brilliant yet simple to use, ageLOC Me delivers five powerful anti-aging products, custom selected from almost 2,000 possible regimen combinations. Delivering one streamlined, simplified regimen to you.
- **ageLOC® LumiSpa® Accent & IdealEyes®**—a LumiSpa attachment and eye cream duo that harnesses the technology of ageLOC LumiSpa to improve skin appearance around the eye area.

FREQUENTLY ASKED QUESTIONS

What is the rationale behind the name?

“Lumi” stands for luminous and emphasizes the initial benefits that using the device and products will have on your skin. “Spa” reinforces the fact that the device, treatment heads, and activating cleansers will help deliver a spa-like beauty experience and results every time you use them.

What are the main differentiators of ageLOC LumiSpa?

There are three key differentiators:

1. Unlike cleansing devices currently on the market, which feature abrasive brush or bristle heads, ageLOC LumiSpa features a gentle, more hygienic silicone treatment head that won't damage the skin or lift skin cells.
2. ageLOC LumiSpa features proprietary Micropulse Oscillation Technology to effectively treat* and cleanse the skin and help deliver superior skin renewal benefits.
3. ageLOC LumiSpa is a dual-action, time-saving device that helps to deliver a superior skin-perfecting beauty experience with deep cleansing benefits.

How does ageLOC LumiSpa work?

With its soft silicone treatment head, ageLOC LumiSpa's proprietary Micropulse Oscillation Technology enables the device to move at the precise frequency needed to promote skin renewal, revealing smoother, softer skin after just one use. At the same time, LumiSpa gently—yet deeply—cleanses away dirt, oil, makeup, and impurities through its proprietary oscillating, pore-purifying action, leaving skin smoother and purer while reducing the appearance of pores.

How do I know which ageLOC LumiSpa treatment head is right for me?

Select your ageLOC LumiSpa treatment head based on your skin type and needs. Each treatment head helps to remove impurities and purify skin, delivering a gentle experience ideal for a variety of skin types. The Normal treatment head is the standard choice and is designed for most skin types. The Gentle treatment head is our gentlest choice and is great for those with sensitive or blemish-prone skin whereas Firm treatment head is for those who want a more intense experience but is not recommended for those with sensitive or blemish-prone skin. You can even personalize and use different heads for different experiences each day.

How long will my treatment head last?

We recommend replacing your treatment head every three months to ensure your skin will receive optimal benefits on a consistent basis.

How does ageLOC LumiSpa work with the activating cleansers?

The activating cleansers were developed specifically for use with ageLOC LumiSpa. They're designed to work synergistically with the oscillating motion of the heads to help optimize distribution of anti-aging ingredients onto the skin while cleansing and exfoliating to help promote healthy, youthful looking skin.

How long should the activating cleanser last?

One activating cleanser should last 30 days.

Are there any patents on ageLOC LumiSpa?

Yes, there are multiple patents issued and pending for ageLOC LumiSpa and the proprietary Micropulse Oscillation Technology.

Do I need to use the ageLOC LumiSpa activating cleansers with ageLOC LumiSpa, or can I use other cleansers with it?

ageLOC LumiSpa activating cleansers are specifically formulated for use with the ageLOC LumiSpa device. They enable the precise cushioning, cleansing, and interaction with the skin, optimizing the effectiveness of LumiSpa's proprietary treatment head action. They are the only formulations tested and approved for use with ageLOC LumiSpa to help deliver the full system benefits and promote healthy, youthful looking skin.

*Treat refers to the delivery of 7 skin benefits: softness, smoothness, radiance, clarity, purified skin, minimized appearance of pores, improved appearance of skin volume & density.
*ageLOC LumiSpa Activating Cleanser-For Blemish Prone is also known as ageLOC LumiSpa Treatment Cleanser-Acne in other markets outside of Southeast Asia.

ageLOC® LumiSpa®

How is ageLOC LumiSpa different from ageLOC Galvanic Spa?

ageLOC LumiSpa provides twice daily treatment* and cleansing and should be used as part of your regular skin care routine, whereas ageLOC Galvanic Spa is a specialized skin care system that can be used three times per week. ageLOC LumiSpa and ageLOC Galvanic Spa deliver complementary anti-aging benefits: ageLOC LumiSpa treats* and cleanses skin for a smoother, softer, more radiant appearance while ageLOC Galvanic Spa conditions, refreshes, and hydrates the skin.

Does ageLOC LumiSpa treat acne?

ageLOC LumiSpa was not designed to treat acne; however, ageLOC LumiSpa Activating Cleanser—For Blemish Prone[^] is a mild formula designed for blemish-prone skin. Because it contains 0.5% salicylic acid, the cleanser helps clear pores and reduce most breakouts.

In a clinical study where subjects used ageLOC LumiSpa with ageLOC LumiSpa Activating Cleanser—For Blemish Prone[^] and Gentle treatment head, statistically significant improvements were noted in as early as one week. These improvements grew to a 71% over the course of a 12-week study. This concludes that ageLOC LumiSpa, ageLOC LumiSpa Activating Cleanser—For Blemish Prone[^] and Gentle treatment head can help mitigate acne breakouts and help the skin look smoother.

We recommend to follow-up with Nu Skin's Clear Action Toner, Day Treatment and Night Treatment to:

- Help improve the appearance of uneven texture and skin discoloration from past breakouts.
- Help reduce the appearance of acne blemishes.
- Help prevent future blemishes from forming.

ageLOC LumiSpa was developed to promote healthy, youthful looking skin and to deliver general skin renewal and purification benefits. Those treating severe acne should continue with their regular skin care regimen and consult a dermatologist before using ageLOC LumiSpa.

Do I need to wash my makeup off before using LumiSpa?

ageLOC LumiSpa conveniently and effectively cleanses away makeup, so you do not need to wash makeup off before using the device. However, to remove heavy makeup, you might want to consider using another cleanser such as Nu Skin Creamy Cleansing Lotion in addition to your LumiSpa. LumiSpa should not be used on the eyes, so

you should use a separate eye makeup remover to remove makeup on and around the eyes.

How do I know if I'm using the right motion and right amount of pressure?

ageLOC LumiSpa has a built-in system that will guide you in using the right motion and right amount of pressure. If you use a scrubbing motion or press too hard, the device will pause and vibrate 3 times until you apply normal speed and pressure. Slowly move the device across your skin in broad strokes. The device will pause briefly every 30 seconds to prompt you to move from one area of the face to the next.

Can I use ageLOC LumiSpa on my body?

ageLOC LumiSpa was designed to target facial skin and deliver anti-aging and deep cleansing benefits in this area. We do not recommend it be used elsewhere on the body.

Can I use ageLOC LumiSpa in the shower?

Yes, ageLOC LumiSpa is waterproof, and you can use it in the shower or other wet environments.

How long does ageLOC LumiSpa last after a full charge?

ageLOC LumiSpa was developed to last at least a week between charges with regular use.

AGELOC LUMISPA ACTIVATING CLEANSERS INGREDIENTS

AGELOC LUMISPA ACTIVATING CLEANSER—NORMAL/COMBO
Water (Aqua), Coco-Betaine, Sodium Cocoyl Glycinate, Glycerin, Acrylates Copolymer, Diisopropyl Sebacate, Glycol Distearate, PEG-8, Rhododendron Ferrugineum Extract, Nannochloropsis Oculata Extract, Schizandra Chinensis Fruit Extract, Porphyrin Umbilicalis Extract, Silybum Marianum Seed Oil, Narcissus Tazetta Bulb Extract, Oryza Sativa (Rice) Extract, Acrylates/C10-30 Alkyl Acrylate Crosspolymer, Polyquaternium-39, Phenoxyethanol, Chlorphenesin, Sodium Benzoate, Phytic Acid, Sodium Hydroxide, Sodium Phytate, Fragrance (Parfum), Synthetic Fluorphlogopite, Titanium Dioxide (CI 77891), Alumina, Tin Oxide.

*Treat refers to the delivery of 7 skin benefits: softness, smoothness, radiance, clarity, purified skin, minimized appearance of pores, improved appearance of skin volume & density.
^ageLOC LumiSpa Activating Cleanser—For Blemish Prone is also known as ageLOC LumiSpa Treatment Cleanser—Acne in other markets outside of Southeast Asia.

**AGELOC LUMISPA ACTIVATING CLEANSER—OILY**

Water (Aqua), Cocamidopropyl Betaine, Acrylates Copolymer, Ammonium Laureth Sulfate, Glycerin, PEG-80 Sorbitan Laurate, Di-PPG-2 Myreth-10 Adipate, Coco-Glucoside, Punica Granatum Juice Extract, Narcissus Tazetta Bulb Extract, Polygonum Cuspidatum Root Extract, Schizandra Chinensis Fruit Extract, Oryza Sativa (Rice) Extract, Glyceryl Oleate, Polyquaternium-10, Myristyl Alcohol, PCA, Carnosine, 1,2-Hexanediol, Fragrance (Parfum), Citric Acid, Aminomethyl Propanol, Dehydroacetic Acid, Chlorphenesin, Phenoxyethanol, Sodium Benzoate, Propanediol, Caprylyl Glycol, Disodium EDTA, Mica (CI 77019), Titanium Dioxide (CI 77891).

AGELOC LUMISPA ACTIVATING CLEANSER—DRY

Water (Aqua), Dicaprylyl Carbonate, Diisostearyl Malate, Tribehenin PEG-20 Esters, Glycerin, Glycereth-26, Squalane, Cetearyl Alcohol, Dimethicone, Caprylic/Capric/Myristic/Stearic Triglyceride, Sodium Astrocaryum Murumuru, Oryza Sativa (Rice) Extract, Schizandra Chinensis Fruit Extract, Xanthophyll, Narcissus Tazetta Bulb Extract, Acrylates/ C10-30 Alkyl Acrylate Crosspolymer, Niacinamide, Urea, Methylsilanol Hydroxyproline Aspartate, PEG-40 Hydrogenated Castor Oil, Sodium Lactate, Sodium PCA, Vegetable Oil, Xanthan Gum, Batyl Alcohol, Butylene Glycol, Pentylene Glycol, Phytosteryl/Octyldodecyl Lauroyl Glutamate, Maltodextrin, Brassicamidopropyl Dimethylamine, Hydroxyphenyl Propamidobenzoic Acid, Salicyloyl Phytosphingosine, Fructose, Glycine, Inositol, Lactic Acid, Lecithin, Phenoxyethanol, Hydroxyacetophenone, Disodium EDTA, Sodium Hydroxide, Fragrance (Parfum).

AGELOC LUMISPA ACTIVATING CLEANSER—SENSITIVE

Water (Aqua), Cetearyl Alcohol, Glycerin, PEG-8 Dimethicone, Cocamidopropyl Betaine, Caprylyl/Capryl Glucoside, Capryl/Capramidopropyl Betaine, Carbomer, Narcissus Tazetta Bulb Extract, Pinus Tabulaeformis Bark Extract, Schizandra Chinensis Fruit Extract, Avena Sativa (Oat) Kernel Extract, Oryza Sativa (Rice) Extract, Allantoin, Zingiber Officinale (Ginger) Root Extract, Bisabolol, Carnosine, Xanthan Gum, Hydroxyacetophenone, Phenoxyethanol, Disodium EDTA, Sodium Hydroxide, Sodium Chloride, Citric Acid.

AGELOC LUMISPA ACTIVATING CLEANSER—FOR BLEMISH PRONE^{*}

Water (Aqua), Cetearyl Alcohol, Glycerin, PEG-8 Dimethicone, Cocamidopropyl Betaine, Caprylyl/Capryl Glucoside, Capryl/Capramidopropyl Betaine, Carbomer, Salicylic Acid, Narcissus Tazetta Bulb Extract, Epilobium Angustifolium Flower/Leaf/Stem Extract, Schizandra Chinensis Fruit Extract, Xanthan Gum, Carnosine, Oryza Sativa (Rice) Extract, Disodium EDTA, Phenoxyethanol, Hydroxyacetophenone, Sodium Hydroxide, Sodium Chloride.

^{*}ageLOC LumiSpa Activating Cleanser-For Blemish Prone is also known as ageLOC LumiSpa Treatment Cleanser-Acne in other markets outside of Southeast Asia.

