




SYSTEM OVERVIEW

Are you ready for a first-of-its-kind, dual-action skin care system that delivers dramatic skin renewal and deep-cleansing in a single treatment? Invented by Nu Skin and dermatologist-tested, ageLOC LumiSpa is inspired by gene expression studies of skin protein production.

With its soft silicone treatment head, ageLOC LumiSpa moves at the precise frequency needed to promote skin renewal, revealing smoother, softer skin after just one use.

LumiSpa gently—yet deeply—cleanses away dirt, oil, makeup, pollutants, and toxins through its proprietary oscillating, pore-tightening action. It leaves skin smoother and purer, while reducing the appearance of pores. And the stimulating, massaging LumiSpa experience will leave you looking forward to your next treatment.

Experience brighter, softer, smoother skin and an energizing facial massage. Simply spend two minutes twice a day to enjoy healthier, youthful looking skin in as little as two weeks.

TARGET AUDIENCE

ageLOC LumiSpa is suitable for any skin-conscious adult seeking ways to minimize the visible signs of aging or to maintain and promote healthy, youthful looking skin.

KEY FEATURES

- **Device**—ageLOC LumiSpa is a waterproof and rechargeable handheld personal skin care device. Used together with the specially formulated LumiSpa treatment cleansers, it simultaneously treats and cleanses the face.
- **Treatment Heads**—ageLOC LumiSpa features two treatment heads: **NORMAL** and **FIRM**. Their gentle, nonabrasive silicone surfaces contain antimicrobial silver and are more hygienic than brushes commonly used with skin care devices. Their sophisticated design helps to maximize product efficacy while making skin more receptive to subsequent products applied to the skin.
- **Treatment Cleansers**—ageLOC LumiSpa features four treatment cleanser options designed for different skin types: **Dry**, **Normal/Combo**, **Oily**, and **Sensitive**. These powerful products are formulated using Nu Skin's proprietary ageLOC ingredient blend. Optimally designed for LumiSpa, they are proven to maximize the skin perfecting and cleansing benefits of the device—revealing healthy, bright, youthful looking skin.

BENEFITS

- Provides instant benefits so skin feels softer and smoother after just one use.
- Delivers seven skin benefits: skin renewal, smoothness, radiance, a more refreshed look,

- purified skin, minimized appearance of pores, and increased appearance of volume and density.
- Benefits intensify over time when consistently used for two minutes twice daily.
- ageLOC LumiSpa uses proprietary Micropulse Oscillation technology, a pore-purifying action that draws out dirt, oil, makeup, pollutants, and toxins and lifts them away without irritating skin while helping to support healthy skin proteins.
- Soft silicone treatment heads are embedded with antimicrobial silver and are more hygienic, less abrasive, and easier to clean than bristle heads commonly found on the market.
- When used in the morning, LumiSpa energizes the skin for a fresher, smoother, more rejuvenated complexion.
- When used at night, LumiSpa deeply cleanses the skin, helping massage and de-stress the skin.
- Removes dead skin cells, promoting dynamic skin renewal and improved skin brightness and smoothness.
- Delivers brighter, healthier, more youthful looking skin in as little as two weeks.
- Gentle enough for twice-daily use.
- Unique sweeping action lifts away loose cells while reducing appearance of pores.
- Goes beyond cleansing to deliver a skin-perfecting treatment through its proprietary oscillating design.
- Doesn't strip skin, but retains skin's natural moisture balance as it cleanses.
- Formulated with Nu Skin's proprietary ageLOC ingredient blend to target the sources of aging and preserve the look of youth.
- ageLOC LumiSpa device is waterproof and can be used in the shower or other wet environment.

CLINICALLY PROVEN BENEFITS

A preliminary 12-week study using ageLOC LumiSpa revealed that individuals who used the

product experienced improvements in skin:

- Smoothness
- Softness
- Texture
- Clarity
- Radiance

Participants noticed several of these improvements after a single use.

USAGE

For optimal results, use ageLOC LumiSpa for two minutes twice daily, as your cleansing step. Morning and night, treat and cleanse with ageLOC LumiSpa treatment cleanser. Then follow with your preferred toner, serum, and moisturizer.

AGELOC® LUMISPA™ TREATMENT CLEANSERS

ageLOC LumiSpa features four treatment cleanser options designed for different skin types: Dry, Normal/Combo, Oily, and Sensitive. These powerful products contain Nu Skin's proprietary ageLOC ingredient blend to target the sources of aging and are formulated to maintain healthy, bright, youthful looking skin.

ageLOC treatment cleansers are specifically formulated to enable the precise cushioning, cleansing, and interaction with the skin, optimizing the effectiveness of LumiSpa's proprietary treatment head action.

**During the Global Preview, only Normal/Combo Cleanser will be available.*

COMPLEMENTARY PRODUCTS

- **ageLOC Galvanic Spa/Nu Skin Facial Spa**—The ageLOC Galvanic Spa uses direct galvanic currents to rejuvenate your complexion and increase skin vibrancy. The Nu Skin Facial Spa uses microcurrents to stimulate and tone the skin.
- **ageLOC Me®**—Creates an entirely new approach to customized anti-aging skin care. Brilliant yet simple to use, ageLOC Me delivers five powerful anti-aging products,

custom selected from almost 2,000 possible regimen combinations. Delivering one streamlined, simplified regimen to you.

- **Tru Face® Line Corrector**—Targets signs of aging with peptides that help soften the appearance of lines around your mouth, eyes, and forehead.
- **ageLOC® Tru Face® Essence Ultra**—Formulated with the power of Ethocyn®* and Nu Skin's revolutionary anti-aging science to target the sources of aging that lead to the loss of firmness.

FREQUENTLY ASKED QUESTIONS

How does ageLOC LumiSpa work?

With its soft silicone treatment head, ageLOC LumiSpa's proprietary Micropulse Oscillation technology enables the device to move at the precise frequency needed to promote skin renewal, revealing smoother, softer skin after just one use. At the same time, LumiSpa gently—yet deeply—cleanses away dirt, oil, makeup, pollutants, and toxins through its proprietary oscillating, pore-tightening action, leaving skin smoother and purer while reducing the appearance of pores.

What are the main differentiators of ageLOC LumiSpa?

There are three key differentiators of ageLOC LumiSpa:

1. Unlike cleansing devices currently on the market, which feature abrasive brush or bristle heads, ageLOC LumiSpa features a gentle, more hygienic silicone treatment head that won't damage the skin or lift skin cells.
2. ageLOC LumiSpa features proprietary Micropulse Oscillation technology to stimulate the skin to more effectively treat and cleanse the skin and deliver superior skin renewal benefits.
3. ageLOC LumiSpa is a dual-action, timesaving device that delivers a superior

skin-perfecting treatment and deep cleansing benefits simultaneously.

What is the rationale behind the name?

"Lumi" stands for luminous and emphasizes the initial benefits that using the device and products will have on your skin. "Spa" reinforces the fact that the device, treatment heads, and treatment cleansers will deliver spa-quality results every time you use them.

How do I know which ageLOC LumiSpa treatment head is right for me?

Select your ageLOC LumiSpa treatment head based on your skin type and treatment needs. Both treatment heads remove impurities and purify skin, delivering a gentle experience ideal for sensitive, dry, oily, and normal to combination skin. The Firm treatment head offers a more aggressive cleansing option.

How long will my treatment head last?

We recommend replacing your treatment head every three months to ensure your skin will receive optimal benefits on a consistent basis.

How does ageLOC LumiSpa work with the treatment cleansers?

The treatment cleansers were developed specifically for use with ageLOC LumiSpa. They're designed to work synergistically with the oscillating motion of the heads to optimize distribution of anti-aging ingredients onto the skin while cleansing and exfoliating to promote healthy, youthful looking skin.

How long should the treatment cleanser last?

One treatment cleanser should last 30 days.

How many studies has Nu Skin conducted for ageLOC LumiSpa?

Nu Skin conducted 11 studies during ageLOC LumiSpa development and will be conducting additional clinical studies to substantiate benefits in the near future.

ageLOC LumiSpa™

Are there any patents on ageLOC LumiSpa?

Yes, there are multiple patents issued and pending for ageLOC LumiSpa and the proprietary Micropulse Oscillation technology.

Do I need to use the ageLOC LumiSpa treatment cleansers with ageLOC LumiSpa, or can I use other treatment formulas or cleansers with it?

ageLOC LumiSpa treatment cleansers are specifically formulated for use with the ageLOC LumiSpa device. They enable the precise cushioning, cleansing, and interaction with the skin, optimizing the effectiveness of LumiSpa's proprietary treatment head action. They are the only formulations tested and approved for use with ageLOC LumiSpa to deliver the full system benefits and promote healthy, youthful looking skin.

How is ageLOC LumiSpa different from ageLOC Galvanic Spa® or Nu Skin® Facial Spa?

ageLOC LumiSpa provides twice daily treatment and cleansing and should be used as part of your regular skin care routine, whereas ageLOC Galvanic Spa and Nu Skin Facial Spa are specialized treatments that can be used three times per week. ageLOC LumiSpa and ageLOC Galvanic Spa/Nu Skin Facial Spa deliver complementary anti-aging benefits: ageLOC LumiSpa treats and cleanses skin for a smoother, softer, more radiant appearance while ageLOC Galvanic Spa conditions, revives, and hydrates the skin, and Nu Skin Facial Spa uses microcurrents to tone and contour the skin.

Does ageLOC LumiSpa treat acne?

ageLOC LumiSpa was not designed to treat acne. It was developed to promote healthy, youthful looking skin and to deliver general skin renewal and purification benefits. Those treating acne should continue with their regular skin care regimen and consult a dermatologist before using ageLOC LumiSpa.

Do I need to wash my makeup off before using LumiSpa?

ageLOC LumiSpa conveniently and effectively cleanses away makeup, so you do not need to wash makeup off before using the device. However, to remove heavy makeup, you might want to consider using another cleanser in addition to your LumiSpa. LumiSpa should not be used on the eyes, so you should use a separate eye makeup remover to remove makeup on and around the eyes.

Can I use ageLOC LumiSpa on my body?

ageLOC LumiSpa was designed to target facial skin and deliver anti-aging and deep cleansing benefits in this area. We do not recommend it be used elsewhere on the body.

Can I use ageLOC LumiSpa in the shower?

Yes, ageLOC LumiSpa is waterproof, and you can use it in the shower or other wet environments.

How long does ageLOC LumiSpa last after a full charge?

ageLOC LumiSpa was developed to last at least a week between charges with regular use.



NU SKIN
DISCOVER THE BEST YOU™



Nu Skin Enterprises Philippines, LLC
15th Floor, Octagon Centre
41 San Miguel Avenue, Pasig City 1605
www.nuskin.com.ph

