

GALVANIC SPA® FACIAL GELS

WITH ageLOC™

Pre-Treat Gel



Treatment Gel



ageLOC™ POSITIONING

Nu Skin's exclusive ageLOC science targets arSuperMarkers—the ultimate sources of aging.

PRODUCT OVERVIEW

Unlock the secret of a younger looking you each time you experience a rejuvenating in-home facial using Galvanic Spa Facial Gels with ageLOC. Transform a tired and stressed complexion into a revitalized, rejuvenated you. Exclusive patented technology gives you an at-home alternative to time consuming and costly spa treatments—in just five minutes a few times a week—with the specially formulated Galvanic Spa® Facial Gels with ageLOC.

Galvanic Spa Facial Gels feature ageLOC—a proprietary science that targets aging at its source. When used with the new ageLOC Edition Galvanic Spa System II, up to five times more ageLOC ingredients are delivered to the skin. The result: skin that immediately feels refreshed, hydrated, and incredibly soft and clean.

An easy-to-do two-step process, the pre-treat step soothes and conditions the skin, binding to impurities and preparing your skin for the treatment step, and the treatment step draws out the impurities targeted in the pre-treat step while delivering ageLOC to the skin.

BENEFITS

Pre-Treat Gel

- Binds to impurities; softens and relaxes pores, allowing them to breathe; dissolves and liquefies excess facial oils.
- Soothes and conditions the skin.

Treatment Gel with ageLOC

- Removes impurities targeted in the pre-treat step.
- ageLOC targets the ultimate sources of aging to preserve the look of youth and reduce the appearance of aging.



- Facilitates recovery from stress; revives and invigorates tired skin by promoting cellular energy.
- Enhances skin hydration, and leaves skin feeling incredibly soft, clean, and refreshed.

Both Gels

- Provide skin rejuvenating results in just five minutes a few times a week that you can see and others will notice.
- Work to restore skin's natural vibrancy by removing impurities and promoting cellular energy.
- Supply beneficial ingredients with delivery enhanced by a galvanic current.

SCIENTIFIC RESEARCH

Nu Skin's continuing commitment to research has led to the discovery of age-related super markers or "arSuperMarkers," components of body chemistry that influence how we age and are the ultimate sources of aging. Nu Skin's exclusive ageLOC science targets these arSuperMarkers. arNOX was the first arSuperMarker Nu Skin classified. Advances in this proprietary research, inspired by leading anti-aging scientists, have led Nu Skin to the discovery of Youth Gene Clusters (YGCs), key arSuperMarkers. YGCs are functional groups of genes whose activity patterns change over time, which can lead to aging. ageLOC science supports resetting these YGCs to more youthful patterns of activity in order to reduce the signs of aging and reveal a more youthful looking you, now and in the future.

KEY INGREDIENTS

Pre-Treat Gel

- Marigold, sea kelp, and vanilla extracts soothe and condition the skin.

GALVANIC SPA® FACIAL GELS WITH ageLOC™

Treatment Gel with ageLOC

- ageLOC compliant blend that targets arSuperMarkers—the ultimate sources of aging.
- Arginine facilitates recovery from stress.
- Magnesium promotes cellular energy to leave skin invigorated and revitalized.

COMPLEMENTARY PRODUCTS

- **ageLOC™ Elements**—the modern essentials of anti-aging skin care, these products purify, moisturize, renew, and reveal younger looking skin.
- **ageLOC™ Transformation**—infuse your skin with the power of ageLOC™ on a daily basis to target aging at its source. These four premier products work together to reveal younger looking skin in eight ways—for a more youthful, healthier looking you now and in the future.
- **ageLOC™ Future Serum**—delivers the maximum concentration of ageLOC to target the ultimate sources of aging. This luxurious serum is clinically proven to reveal younger looking skin in eight significant ways: by promoting youthful skin structure, smoothing texture, reducing the appearance of fine lines and wrinkles, increasing radiance, minimizing the appearance of pores, evening skin tone, reducing discoloration, and increasing hydration.
- **Tru Face® Essence Ultra**—featuring 25 percent more serum than original Tru Face Essence, and combined for the first time with the coenzyme Q10, Tēgreen 97,® and a protective antioxidant network, the powerful science of Tru Face Essence Ultra helps restore facial firmness, promotes cellular energy production, and helps prevent future signs of aging.
- **Tru Face® Line Corrector**—features a revolutionary pro-collagen peptide that helps soften moderate to deep lines around the mouth, eyes, and forehead.
- **Nu Skin 180°® Anti-Aging Skin Therapy System**—designed to reverse the full range of aging effects, each product in this system plays an integral role in delivering anti-aging benefits and works synergistically to enhance results.

FREQUENTLY ASKED QUESTIONS

What is new about the Galvanic Spa Facial Gels with ageLOC?
Nu Skin's ageLOC science targets the ultimate sources of aging, arSuperMarkers. The proprietary ingredients contained in Galvanic Spa Facial Gels with ageLOC have been validated with new genomic research showing that ageLOC ingredients help reset Youth Gene Clusters to reflect a more youthful gene expression profile in skin.

Can people with medical conditions use the Galvanic Spa System II?

Check with your physician before using the Galvanic Spa II instrument if you are pregnant, have a pacemaker or similar device, are epileptic, wear metal braces on your teeth, or have a metal implant. Also, be sure to check the instructions and directions for use.

Can people with acne use Galvanic Spa Facial Gels with ageLOC?

Galvanic Spa Facial Gels with ageLOC are great for any skin type. They should not cause any additional problems for individuals with acne, but are not designed to help or treat this condition. Additionally, the Galvanic Spa II instrument should not be used over an area with many broken capillaries or over an inflamed or open wound.

Can I use part of one gel for one application and then use the rest a few days later?

The facial gels are packaged in individually metered dosages. Using a full vial as directed will ensure that each application provides the most beneficial amount of key ingredients and allows for optimal conductivity for the galvanic current. We recommend that the facial gels be used in one sitting and not be saved for later use.

How often should I use the facial gels?

You can use the facial gels up to two or three times a week.

Can I use another kind of gel, like Enhancer, with the Galvanic Spa II instrument in place of the facial gels?

Using a non-charged gel, such as Enhancer, with the Galvanic Spa would not provide any added benefit. The Galvanic Spa is designed for use with the Galvanic Spa Facial Gels with ageLOC and other Galvanic Spa compatible products, which are specifically formulated to work with the galvanic current and have been tested for safety with the Galvanic Spa instrument. Other products not specifically formulated to work with the Galvanic Spa may cause irritation.

INGREDIENTS

Galvanic Spa Pre-Treat Gel

Water (Aqua), Butylene Glycol, Glycerin, Calendula Officinalis Flower Extract, Macrocyctis Pyrifera Extract, Vanilla Planifolia Fruit Extract, Hydroxypropyl Methylcellulose, Glyceryl Acrylate/Acrylic Acid Copolymer, Polysorbate 20, TEA-Cocoyl Glutamate, Propylene Glycol,

GALVANIC SPA® FACIAL GELS WITH ageLOC™

Carbomer, Potassium Hydroxide, Citric Acid, Disodium EDTA, Chlorphenesin, Phenoxyethanol.

Galvanic Spa Treatment Gel

Water (Aqua), Butylene Glycol, Glycerin, Isoceteth-20, Schizandra Chinensis Fruit Extract, Narcissus Tazetta Bulb Extract, Tocopheryl Nicotinate, Arginine, Magnesium Aspartate, Salicin, Jasminum Officinale (Jasmine) Extract, Rosa Damascena Flower Oil, Santalum Album (Sandalwood) Oil, TEA-Cocoyl Glutamate, Polysorbate 20, Oleth-10, Hydroxypropyl Methylcellulose, Xanthan Gum, Glyceryl Acrylate/Acrylic Acid Copolymer, Propylene Glycol, Disodium EDTA, Citric Acid, Chlorphenesin, Potassium Sorbate, Phenoxyethanol.

