

ageLOC® LumiSpa®



SYSTEM OVERVIEW

Are you ready for a first-of-its-kind, dual-action skin care system that delivers dramatic skin renewal and deep-cleansing in a single treatment? Developed in conjunction with a board-certified dermatologist, ageLOC LumiSpa is inspired by gene expression studies of skin protein production.

With its gentle and soft silicone treatment head, ageLOC LumiSpa moves at the precise frequency to promote skin renewal, revealing smoother, softer skin after just one use.

ageLOC LumiSpa gently—yet deeply—cleanses away dirt, oil, makeup, and toxins through its patented counter-rotating, pore-tightening action. It leaves skin smoother and purer, while reducing the appearance of pores. And the stimulating, massaging ageLOC LumiSpa experience will leave you looking forward to your next treatment.

Experience brighter, softer, smoother skin and an energizing facial massage. Simply spend two minutes twice a day to enjoy healthier, youthful looking skin in as little as two weeks.

TARGET AUDIENCE

ageLOC LumiSpa is suitable for any skin-conscious adult seeking ways to minimize the visible signs of aging or to maintain and promote healthy, youthful looking skin.

KEY FEATURES

- **Device**—ageLOC LumiSpa is a waterproof and rechargeable handheld personal skin care device. Used together with the specially formulated ageLOC LumiSpa Activating Cleansers, it simultaneously treats and cleanses the face.

- **Treatment Heads**—ageLOC LumiSpa features two treatment heads: NORMAL and FIRM. Their gentle, non abrasive silicone surfaces contain antimicrobial silver and is more hygienic than brushes commonly used with other devices. Their sophisticated design helps to maximize product efficacy while making skin more receptive to subsequent products applied to the skin.
- **Activating Cleansers**—ageLOC LumiSpa features four Activating Cleanser options designed for different skin types:
 - Dry skin
 - Normal/Combo skin
 - Oily skin
 - Sensitive skin
 These powerful products are formulated using Nu Skin's proprietary ageLOC ingredient blend. Optimally designed for ageLOC LumiSpa, they are proven to maximize the skin perfecting and cleansing benefits of the device—revealing healthy, bright, youthful looking skin.

BENEFITS

- Provides instant benefits so skin feels softer and smoother after just one use.
- Delivers seven skin benefits:
 - Skin renewal
 - Smoothness
 - Radiance
 - A more refreshed look
 - Purified skin
 - Minimized appearance of pores
 - Increased appearance of volume and density
- ageLOC LumiSpa benefits intensify over time when consistently used for two minutes twice daily.
- When used in the morning, ageLOC LumiSpa energizes the skin for a fresher, smoother, more rejuvenated complexion.

- When used in the evening, ageLOC LumiSpa deeply cleanses the skin, and helps to massage and de-stress the skin.
- Removes dead skin cells, promoting dynamic skin renewal and improved skin brightness and smoothness.
- Soft silicone treatment heads are more hygienic, less abrasive, and easier to clean than bristle heads commonly found in the market.
- Delivers brighter, healthier, more youthful looking skin in as little as two weeks.
- Gentle enough for twice-daily use.
- Unique sweeping action lifts away loose cells while reducing the appearance of pores.
- Goes beyond cleansing to deliver a skin-perfecting treatment through its proprietary oscillation design.
- Gentle and does not strip skin, but retains skin's natural moisture balance as it cleanses.
- Gently cleanses away impurities, makeup, dirt, and toxins without irritating skin.
- Formulated with Nu Skin's proprietary ageLOC ingredient blend to target the sources of aging and preserve the look of youth.

CLINICALLY PROVEN BENEFITS

A preliminary 12-week study using ageLOC LumiSpa revealed that individuals who used the product experienced improvements in skin:

- Smoothness
- Softness
- Texture
- Clarity
- Radiance

Participants noticed several of these improvements after a single use.

USAGE

For optimal results, use ageLOC LumiSpa for two minutes twice daily as your cleansing step. Morning and evening, treat and cleanse with ageLOC LumiSpa Activating Cleanser. Then followed by your preferred toner, serum, and moisturizer.

ageLOC® LUMISPA® ACTIVATING CLEANSERS

ageLOC LumiSpa features four Activating Cleansers options designed for different skin types: dry, normal/ combo, oily, and sensitive skin. These powerful products contain Nu Skin's proprietary ageLOC ingredient blend to target the sources of aging and are formulated to maintain healthy, bright, youthful looking skin.

ageLOC LumiSpa Activating Cleansers are specifically formulated to enable precise cushioning, cleansing, and interaction with the skin, optimizing the effectiveness of ageLOC LumiSpa's proprietary treatment head action.

ageLOC® LUMISPA® ACTIVATING CLEANSERS KEY INGREDIENTS

- **ageLOC LumiSpa Activating Cleanser—Dry**
 - Contains Sodium Lactate, Sodium PCA, Glycine, Fructose, Niacinamide, and Inositol to help protect skin's natural moisture.
 - Squalane to help maintain skin's natural moisture barrier.
 - Salicyloyl Phytosphingosine to soothe the skin.
- **ageLOC LumiSpa Activating Cleanser—Normal/ Combo**
 - Contains Milk Thistle (Silybum Marianum Seed) Oil, Porphyra Umbilicalis Extract, and Nannochloropsis Oculata Extract to help maintain skin's natural moisture barrier.
 - Rhododendron (Ferrugineum) Extract to provide antioxidant protection.
- **ageLOC LumiSpa Activating Cleanser—Oily**
 - Contains Japanese Knotweed (Polygonum Cuspidatum Root) Extract and Myristyl PCA to help reduce oil production and soothe the skin.
 - Pomegranate (Punica Granatum) Juice Extract and Carnosine to provide antioxidant protection.
- **ageLOC LumiSpa Activating Cleanser—Sensitive**
 - Most gentle formula. It is a pink gel cream formulated with a mild surfactant system, which results in a low foam.
 - It contains Avena Sativa (Oat) Kernel Extract, Pinus Tabulaeformis Bark Extract (which contributes to the pink color) and Bisabolol to help calm and soothe the delicate skin.
 - Allantoin helps to condition the skin.

COMPLEMENTARY PRODUCTS

- **ageLOC® Galvanic Spa®**—The ageLOC Galvanic Spa utilizes direct/galvanic currents to rejuvenate your complexion and increase skin vibrancy. Use ageLOC Galvanic Spa two to three times weekly along with daily ageLOC LumiSpa treatments to dramatically increase visible results.
- **ageLOC Me®**—Creates an entirely new approach to customized anti-aging skin care. Brilliant yet simple to use, ageLOC Me delivers five powerful anti-aging products, custom selected from more than 2,000 possible regimen combinations. Delivering one streamlined, simplified regimen to you.
- **Tru Face® Line Corrector**—Targets signs of aging with pro-collagen peptides that helps soften moderate to deep lines around your mouth, eyes, and forehead. Scientists understand that peptides have the ability to

ageLOC® LumiSpa®

send age-reversing signals to collagen-producing cells, promoting younger looking skin.

- **ageLOC® Tru Face® Essence Ultra**—Formulated with the power of Ethocyn®* and Nu Skin's revolutionary anti-aging science to target the sources of aging that lead to the loss of firmness. Leave undefined skin in the past and enjoy firmer, younger looking skin today and in the future with ageLOC Tru Face Essence Ultra.

FREQUENTLY ASKED QUESTIONS

How does ageLOC LumiSpa work?

With its soft silicone treatment head, ageLOC LumiSpa's proprietary Micropulse Oscillation technology enables the device to move at the precise frequency to promote skin renewal, revealing smoother, softer skin after just one use. At the same time, ageLOC LumiSpa gently—yet deeply—cleanses away dirt, oil, makeup, and toxins through its proprietary oscillation, pore-tightening action, leaving skin smoother and purer while reducing the appearance of pores.

What are the main differentiators of ageLOC LumiSpa?

There are three key differentiators of ageLOC LumiSpa:

1. Unlike cleansing devices currently in the market, which feature abrasive brush or bristle heads, ageLOC LumiSpa features a gentle, more hygienic silicone treatment head that won't damage the skin or lift skin cells.
2. ageLOC LumiSpa features patented Micropulse Oscillation technology to stimulate the skin to more effectively treat and cleanse the skin and deliver superior skin renewal benefits.
3. ageLOC LumiSpa is a dual-action, time-saving device that delivers a superior skin-perfecting treatment and deep cleansing benefits simultaneously.

Do I need to remove my makeup before using ageLOC LumiSpa?

ageLOC LumiSpa conveniently and effectively cleanses away makeup, so you do not need to remove makeup before using LumiSpa. However, to remove heavy makeup, you might want to consider using a makeup remover before your LumiSpa. LumiSpa should not be used on the eyes, so you should use a separate eye makeup remover to remove makeup on and around the eyes.

How do I know which ageLOC LumiSpa treatment head is right for me?

You can select your ageLOC LumiSpa treatment head based on your skin type and treatment needs. The Normal treatment head removes impurities and purifies skin, delivering a gentle experience ideal for sensitive, dry, and normal to combination skin. The Firm treatment head exfoliates skin for a deeper cleanse ideal for oily skin.

How long will my treatment head last?

We recommend replacing your treatment head every three months to ensure your skin will receive optimal benefits on a consistent basis.

How does ageLOC LumiSpa work with the Activating Cleansers?

The Activating Cleansers are developed specifically for use with ageLOC LumiSpa. They are designed to work synergistically with the counter-rotating motion of the heads to optimize distribution of anti-aging ingredients onto the skin while cleansing and exfoliating to promote healthy, youthful looking skin.

How long can the ageLOC LumiSpa Activating cleanser last?

Each Activating Cleanser should last 30 days, based on twice daily usage.

How many studies has Nu Skin conducted for ageLOC LumiSpa?

Nu Skin conducted 11 studies during ageLOC LumiSpa development and will be conducting additional clinical studies to substantiate benefits in the near future.

Are there any patents on ageLOC LumiSpa?

Yes, there are multiple patents issued and pending for ageLOC LumiSpa and the proprietary Micropulse Oscillation technology.

Do I need to use the ageLOC LumiSpa Activating Cleansers with ageLOC LumiSpa, or can I use other treatment formulas or cleansers with it?

ageLOC LumiSpa Activating Cleansers are specifically formulated for use with the ageLOC LumiSpa device. They enable the precise cushioning, cleansing, and interaction with the skin, optimizing the effectiveness of LumiSpa's proprietary treatment head action. They are the only formulations tested and approved for use with ageLOC LumiSpa to deliver the full system benefits and promote healthy, youthful looking skin.

How is ageLOC LumiSpa different from ageLOC Galvanic Spa?

ageLOC LumiSpa provides twice daily treatment and cleansing benefits and should be used as part of your regular skin care routine, whereas ageLOC Galvanic Spa is a specialized treatment that can be used up to three times per week. ageLOC LumiSpa and ageLOC Galvanic Spa delivers complementary anti-aging benefits: ageLOC LumiSpa treats and cleanses skin for a smoother, softer, more radiant appearance while ageLOC Galvanic Spa conditions, revives, and hydrates the skin.

ageLOC® LumiSpa®

Can ageLOC LumiSpa treat acne?

ageLOC LumiSpa is not designed to treat acne. The system is developed to promote healthy, youthful looking skin and to deliver general skin renewal and purification benefits. Those treating acne should continue with their regular skin care regimen and consult a dermatologist before using ageLOC LumiSpa.

Can I use ageLOC LumiSpa on my body?

ageLOC LumiSpa is designed to target facial skin and deliver anti-aging and deep cleansing benefits in this area. We do not recommend it to be used on the body.

Can I use ageLOC LumiSpa in the shower?

Yes, ageLOC LumiSpa is waterproof, so you can use it in the shower or other wet environments. Please note that the charging base is not waterproof and it should not come into contact with water.

How long does ageLOC LumiSpa last after a full charge?

ageLOC LumiSpa is developed to last for approximately a week between charges with regular use.

INGREDIENTS

ageLOC LUMISPA ACTIVATING CLEANSER—DRY

Water (Aqua), Dicaprylyl Carbonate, Diisostearyl Malate, Tribehenin PEG-20 Esters, Glycerin, Glycereth-26, Squalane, Cetearyl Alcohol, Dimethicone, Caprylic/Capric/Myristic/Stearic Triglyceride, Sodium Astrocaryum Murumurate, Phenoxyethanol, Butylene Glycol, Pentylene Glycol, Acrylates/C10-30 Alkyl Acrylate Crosspolymer, PEG-40 Hydrogenated Castor Oil, Sodium PCA, Sodium Lactate, Hydroxyacetophenone, Vegetable Oil, Xanthan Gum, Batyl Alcohol, Fragrance (Parfum), Methylsilanol Hydroxyproline Aspartate, Phytosteryl/Octyldodecyl Lauroyl Glutamate, Maltodextrin, Sodium Hydroxide, Salicyloyl Phytosphingosine, Disodium EDTA, Brassicamidopropyl Dimethylamine, Hydroxyphenyl Propamidobenzoic Acid, Oryza Sativa (Rice) Extract, Schizandra Chinensis Fruit Extract, Sodium Benzoate, Inositol, Urea, Glycine, Lactic Acid, Fructose, Niacinamide, Narcissus Tazetta Bulb Extract, Xanthophyll, Lecithin.

ageLOC LUMISPA ACTIVATING CLEANSER—NORMAL/COMBO

Water (Aqua), Coco-Betaine, Sodium Cocoyl Glycinate, Glycerin, Acrylates Copolymer, Diisopropyl Sebacate, Glycol Distearate, Phenoxyethanol, PEG-8, Acrylates/C10-30 Alkyl Acrylate Crosspolymer, Polyquaternium-39, Sodium Benzoate, Chlorphenesin, Fragrance (Parfum), Synthetic Fluorphlogopite, Silybum Marianum Seed Oil,

Sodium Hydroxide, Sodium Phytate, Titanium Dioxide (CI 77891), Phytic Acid, Oryza Sativa (Rice) Extract, Schizandra Chinensis Fruit Extract, Alumina, Tin Oxide, Nannochloropsis Oculata Extract, Porphyra Umbilicalis Extract, Narcissus Tazetta Bulb Extract, Rhododendron Ferrugineum Extract.

ageLOC LUMISPA ACTIVATING CLEANSER—OILY

Water (Aqua), Cocamidopropyl Betaine, Acrylates Copolymer, Ammonium Laureth Sulfate, Glycerin, PEG-80 Sorbitan Laurate, Di-PPG-2 Myreth-10 Adipate, Aminomethyl Propanol, Glyceryl Oleate, Coco-Glucoside, Phenoxyethanol, 1,2-Hexanediol, Caprylyl Glycol, Chlorphenesin, Sodium Benzoate, Disodium EDTA, Polyquaternium-10, Dehydroacetic Acid, Fragrance (Parfum), Titanium Dioxide (CI 77891), Mica (CI 77019), Propanediol, Oryza Sativa (Rice) Extract, Myristyl Alcohol, PCA, Schizandra Chinensis Fruit Extract, Carnosine, Punica Granatum Juice Extract, Narcissus Tazetta Bulb Extract, Polygonum Cuspidatum Root Extract, Citric Acid.

ageLOC LUMISPA ACTIVATING CLEANSER—SENSITIVE

Water (Aqua), Cetearyl Alcohol, Glycerin, PEG-8 Dimethicone, Cocamidopropyl Betaine, Caprylyl/Capryl Glucoside, Capryl/Capramidopropyl Betaine, Carbomer, Phenoxyethanol, Hydroxyacetophenone, Xanthan Gum, Sodium Chloride, Bisabolol, Carnosine, Disodium EDTA, Oryza Sativa (Rice) Extract, Pinus Tabulaeformis Bark Extract, Allantoin, Avena Sativa (Oat) Kernel Extract, Schizandra Chinensis Fruit Extract, Narcissus Tazetta Bulb Extract, Zingiber Ocinalle (Ginger) Root Extract, Sodium Hydroxide, Citric Acid.